



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Nutrient Agar w/ 1% Peptone

Product Code: DM 1012

Application: Nutrient agar medium is used as a general purpose culture media which may be used as enriched medium by incorporating 10% blood or other biological fluid.

Composition**

Ingredients	Gms / Litre
Peptic digest of animal tissue	10.000
Sodium chloride	5.000
Beef extract	5.000
Agar	15.000
Final pH (at 25°C)	7.4±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Nutrient Agar with 1% Peptone has almost double concentration of the nitrogen sources than that used in Nutrient Agar, making media more nutritive. Nutrient media are basic culture media used for maintaining microorganisms, cultivating fastidious organisms by enriching media with serum or blood and are also used for purity checking prior to biochemical or serological testing^(1, 2). It is one of the several non-selective media useful in routine cultivation of microorganisms^(3, 4).

Beef extract and peptic digest of animal tissue provide the necessary nitrogen compounds, carbon, vitamins and also some trace ingredients to non-fastidious organisms like *Bacillus subtilis* and *Staphylococcus aureus*. Sodium chloride maintains osmotic equilibrium of the medium. With the addition of 10% v/v blood or other biological fluids like ascetic fluid, serum etc, this media is recommended for growing fastidious organisms.

Peptic digest of animal tissue and beef extract are nutritionally rich in supplying essential nitrogen and growth factors⁽⁵⁾.

Methodology

Suspend 35 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well before pouring.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Basal medium: Light yellow coloured clear to slightly opalescent gel. After addition of 5%v/v sterile defibrinated blood : Cherry red coloured opaque gel forms in Petri plates

Reaction

Reaction of 3.5% w/v aqueous solution at 25°C. pH : 7.4±0.2

pH Range:-

7.20-7.60

Cultural Response/Characteristics

DM 1012: Cultural characteristics observed with added 5% v/v sterile defibrinated blood, after an incubation at 35-37°C for 18-48 hours.





Dehydrated Culture Media
Bases / Media Supplements

Organism	Inoculum (CFU)	Growth	Recovery	Growth w/ Blood	Recovery w/ Blood	Haemolysis
<i>Neisseria meningitidis</i> ATCC 13090	50-100	good	50-70%	luxuriant	>=70%	none
<i>Staphylococcus aureus</i> ATCC 25923	50-100	luxuriant	>=70%	luxuriant	>=70%	beta
<i>Streptococcus pneumoniae</i> ATCC 6303	50-100	good	50-70%	luxuriant	>=70%	alpha
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	good	50-70%	luxuriant	>=70%	beta
<i>Staphylococcus aureus</i> ATCC 6538	50-100	luxuriant	>=70%	luxuriant	>=70%	beta

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Lapage S., Shelton J. and Mitchell T., 1970, Methods in Microbiology', Norris J. and Ribbons D., (Eds.), Vol. 3A, Academic Press, London.
2. MacFaddin J. F., 2000, Biochemical Tests for Identification of Medical Bacteria, 3rd Ed., Lippincott, Williams and Wilkins, Baltimore.
3. Downes F. P. and Ito K., (Ed.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., American Public Health Association, Washington, D.C.
4. American Public Health Association, Standard Methods for the Examination of Dairy Products, 1978, 14th Ed., Washington D.C.
5. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. 1, Williams and Wilkins, Baltimore

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specificatons for identity and performens parameters.

