

Technical Information

Phenol Red Broth Base

Product Code: DM 1054

Application: - Phenol Red Broth Base is used as a basal medium to which carbohydrates are added for determination of fermentation reactions of pure cultures of microorganisms.

Composition**

Ingredients	Gms / Litre
Proteose peptone	10.000
Beef extract	1.000
Sodium chloride	5.000
Phenol red	0.018
Final pH (at 25°C)	7.4±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Phenol Red Broth Medium devised by Vera ⁽²⁾ is recommended to determine the fermentation reaction of carbohydrates for the differentiation of microorganisms ⁽³⁻⁵⁾. Phenol Red Broth Medium with different of added carbohydrates serves as a differential medium by helping in differentiation of different species and genera by their ability to ferment the specific carbohydrate, with the production of acid or acid and gas ⁽⁶⁾.

Phenol Red Broth Base is a complete medium without added carbohydrate, which can be used with the addition of 5-10 %, desired carbohydrate. It is used as a negative control for studying fermentations or as a base for the addition of carbohydrates. Proteose peptone and beef extract serve as sources for carbon and nitrogen. Sodium chloride is the osmotic stabilizer. Phenol red is the pH indicator, which turns yellow at acidic pH. Gas formation is seen in Durhams tubes. All of the member of family *Enterobacteriaceae* grows well in this medium. In addition to producing a pH colour shift, the production of mixed acids, notably butyric acids, often results in a pungent, foul odour from the culture medium ⁽¹⁾.

Methodology

Suspend 16 grams of powder media in 1000 ml distilled water, mix well. Shake well & heat if necessary to ensure complete solution. Distribute in fermentation tubes (tubes containing inverted Durham's tubes). Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Aseptically add filter sterilized or autoclave sterilized carbohydrate solution to sterile basal medium.

Quality Control

Physical Appearance

Light yellow to pink coloured homogeneous free flowing powder

Colour and Clarity of prepared medium

Red coloured clear solution without any precipitate

Reaction

Reaction of 1.6% w/v aqueous solutions at 25°C. pH : 7.4±0.2

pH Range 7.20-7.60

Cultural Response/Characteristics

DM 1054: Cultural characteristics observed after an incubation at 35 - 37°C for 18 - 24 hours.(longer if necessary)

Organism	Inoculum (CFU)	Growth	Without carbohydrate (Acid)	Without carbohydrate (Gas)	With dextrose (Acid)	With dextrose (Gas)
----------	----------------	--------	-----------------------------	----------------------------	----------------------	---------------------

<i>Citrobacter freundii</i> ATCC 8090	50-100	luxuriant	Negative, reaction no colour change	Negative reaction	Positive reaction, yellow colour	Positive reaction
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	Negative reaction, no colour change	Negative reaction	Positive reaction, yellow colour	Positive reaction
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	Negative reaction, no colour change	Negative reaction	Positive reaction, yellow colour	Positive reaction
<i>Klebsiella pneumoniae</i> ATCC 13883	50-100	luxuriant	Negative reaction, no colour change	Negative reaction	Positive reaction, yellow colour	Positive reaction
<i>Proteus vulgaris</i> ATCC 13315	50-100	luxuriant	Negative reaction, no colour change	Negative reaction	Positive reaction, yellow colour	Positive reaction
<i>Salmonella Typhi</i> ATCC 6539	50-100	luxuriant	Negative reaction, no colour change	Negative reaction	Positive reaction, yellow colour	Positive reaction
<i>Salmonella Typhimurium</i> ATCC 14028	50-100	luxuriant	Negative reaction, no colour change	Negative reaction	Positive reaction, yellow colour	Positive reaction
<i>Serratia marcescens</i> ATCC 8100	50-100	luxuriant	Negative reaction, no colour change	Negative reaction	Positive reaction, yellow colour	Positive reaction
<i>Shigella flexneri</i> ATCC 12022	50-100	luxuriant	Negative reaction, no colour change	Negative reaction	Positive reaction, yellow colour	Negative reaction

Storage and Shelf Life

Dried media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8⁰ in sealable plastic bags for 2-5 days.

Further Reading

1. Thom and Church, 1926, The Aspergilli, 39.
2. Eaton A. D., Clesceri L. S. and Greenberg A. E., (Ed.), 1998, Standard Methods for the Examination of Water and Waste water, 20th Ed., American Public Health Association. Washington, D.C

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.