



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

MacConkey Agar w/ Bromo Thymol Blue

Product Code: DM 1061

Application: - MacConkey Agar w/ Bromo Thymol Blue is recommended for detection of lactose fermenting enteric bacteria.

Composition**

Ingredients	Gms / Litre
Peptic digest of animal tissue	17.000
Proteose peptone	3.000
Lactose	10.000
Bile salts	1.500
Sodium chloride	5.000
Bromothymol blue	0.030
Agar	15.000
Final pH (25°C)	7.1±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

MacConkey Agar is one of the earliest selective and differential medium for cultivation of enteric microorganisms from a different type of clinical specimens^(1, 2). The original medium contains protein, bile salts, sodium chloride and two dyes. MacConkey Agar w/ Bromo thymol blue is a modification of the original medium by replacing two dyes with a single dye i.e. bromo thymol blue.

Peptic digest of animal tissue and proteose peptone serve as the sources of essential nutrients. Lactose is the fermentable carbohydrate source. Lactose fermenting enteric bacteria ferment lactose and produce acidic byproducts. This acidic condition formed make the pH indicator dye i.e. bromo thymol blue to change colour from blue to yellow. Lactose non-fermenters fail to make a colour change in the medium. Sodium chloride maintains the osmotic equilibrium of the medium. Bile salts serves to make the medium selective by inhibiting the accompanying gram-positive bacteria.

Methodology

Suspend 51.53 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure for 15 minutes. Mix well and pour into sterile Petri plates.

Quality Control

Physical Appearance

Cream to greenish yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Green coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 5.15% w/v aqueous solution at 25°C, pH:-7.1±0.2

pH range 6.90-7.30

Cultural Response/ characteristics

DM 1061: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.





Dehydrated Culture Media
Bases / Media Supplements

Organism	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	>=50 %	yellow
<i>Enterococcus faecalis</i> ATCC 29212	>=10 ³	inhibited	0 %	
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	>=50 %	yellow
<i>Proteus vulgaris</i> ATCC 13315	50-100	luxuriant	>=50 %	colourless-light blue
<i>Salmonella Typhi</i> ATCC 6539	50-100	luxuriant	>=50 %	colourless-light blue
<i>Shigella flexneri</i> ATCC 12022	50-100	luxuriant	>=50 %	colourless-light blue
<i>Staphylococcus aureus</i> ATCC 25923	>=10 ³	inhibited	0 %	

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. MacConkey, 1900, The Lancet, ii:20.
2. MacConkey, 1905, J. Hyg., 5:333.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.

