

## Technical Information

### Deoxycholate Lactose Agar

**Product Code: DM 1066**

**Application:** Deoxycholate Lactose Agar is a differential and slightly selective medium used for the isolation and enumeration of coliforms in water, wastewater, milk and dairy products.

#### Composition\*\*

Ingredients	Gms / Litre
Peptone, special	10.000
Lactose	10.000
Sodium chloride	5.000
Sodium citrate	2.000
Sodium deoxycholate	0.500
Neutral red	0.030
Agar	15.000
Final pH ( at 25°C)	7.1±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

Deoxycholate Lactose Agar is a modification of Deoxycholate Agar as described by Leifson <sup>(1)</sup> and prepared according to formula specified in Standard Methods for Examination of Dairy Products <sup>(2)</sup> Water and Waste Water <sup>(3)</sup> and Food <sup>(4)</sup> for the detection of coliform bacilli. It differs from Deoxycholate Agar (DM1030) due to decreased concentration of sodium deoxycholate. Pour plate method is carried out using suitable dilutions. A thin layer of additional agar can be poured over the solidified pour plates to facilitate enumeration. Deoxycholate Lactose Agar is selective against gram-positive organisms which are inhibited by optimum concentration of sodium deoxycholate and sodium citrate in the medium. It helps to differentiate between lactose fermenting and nonfermenting enteric bacilli. Lactose helps in differentiating enteric bacilli, as lactose fermenters produce red colonies while lactose nonfermenters produce colourless colonies. Coliform bacteria, if present form pink colonies on this medium. The degradation of lactose causes acidification of the medium surrounding the relevant colonies and the pH indicator neutral red changes its colour to red. These colonies usually are also surrounded by a turbid zone of precipitated deoxycholic acid due to acidification of the medium. Sodium deoxycholate combines with neutral red in an acidic environment, causing the dye to go out of the solution with the subsequent precipitation of deoxycholate <sup>(1)</sup>.

#### Methodology

Suspend 42.53 grams of powder media in 1000 ml distilled water. Mix well and heat to dissolve the medium completely. The medium requires no autoclaving if it is to be used at once. If the medium is to be stored, it should be sterilized at 15 lbs pressure (121°C) for 15 minutes. AVOID OVERHEATING.

## Quality Control

### Physical Appearance

Light yellow to pink homogeneous free flowing powder

### Gelling

Firm, comparable with 1.5% Agar gel

### Colour and Clarity of prepared medium

Reddish orange coloured, clear to slightly opalescent gel forms in Petri plates

### Reaction

Reaction of 4.25% w/v aqueous solution at 25°C. pH : 7.1±0.2

**pH Range** 6.90-7.30

### Cultural Response/Characteristics

DM1066: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Colour of Colony
<i>Bacillus subtilis</i> ATCC 6633	$\geq 10^3$	Inhibited	0%	
<i>Escherichia coli</i> ATCC 25922		good-luxuriant	$\geq 50\%$	pink w/bile precipitate
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	good-luxuriant	$\geq 50\%$	Pink
<i>Enterococcus faecalis</i> ATCC 29212	$\geq 10^3$	Inhibited	0%	
<i>Salmonella Typhimurium</i> ATCC 14028	50-100	good-luxuriant	$\geq 50\%$	colourless

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. Leifson, 1935, J. Path. Bact., 40:58.

2. Richardson (Ed.), 1985, Standard Methods for the Examination of Dairy Products, 15th ed., APHA, Washington, D.C. 3. Greenberg A. E., Eaton A. D., Clesceri L. S., (Eds.), 1998, Standard Methods for the Examination of Water and Waste Water, 20th Ed., APHA, Washington, D.C.

4. Speck M. (Ed.), 1984, Compendium of Methods for the Microbiological Examination of Foods, 2nd ed., APHA, Washington, D.C.

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.