

Technical Information

Bile Broth Base

Product Code: DM 1071

Application: Bile Broth Base is used for the cultivation of members of *Enterobacteriaceae* and in blood clots culture from patients with suspected enteric fever.

Composition**

Ingredients	Gms / Litre
Peptic digest of animal tissue	20.000
Sodium taurocholate	5.000
Sodium chloride	5.000
Final pH (at 25°C)	7.6±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

When blood samples from a patient with suspected enteric fever is submitted for the widal test, it is useful as a routine to do the clot culture after separation of serum ⁽¹⁾. If the blood has been withdrawn with strict aseptic precautions, the clot may be placed in a wide tube half-filled with broth, or in a wide mouth screw-capped bottle containing 80 ml of broth. When there is any doubt regarding the presence of contaminating organisms, and there is always a possibility when blood specimens are sent to the laboratory from a distance, the clot should be transferred directly to a tube of sterile ox bile broth and disintegrated with aseptic precautions. After overnight incubation the bile culture is examined for enteric organism in the usual manner. A simple & easy method of clot culture with Streptokinase has been recommended. In this procedure Blood is allowed to clot in 5 ml quantities in sterile screw-capped universal containers. The separated serum is removed and 15 ml of 0.5% Bile Broth Base with Streptokinase 100 units/ml is added to each bottle. The streptokinase causes rapid clot lysis with release of bacteria trapped in the clot ⁽²⁾ Peptic digest of animal tissue serves important sources of nitrogen. Sodium taurocholate inhibits majority of Gram-positive species. Sodium chloride maintains the isotonicity of the medium whereas addition of streptokinase solution causes rapid clot lysis with release of bacteria trapped in the clot.

Methodology

Suspend 30.0 grams of powder media in 1000 ml distilled water. Shake well & heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 40°C and add 1 ml of Streptokinase solution (100000 units/ml). Mix well and dispense as desired.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Yellow coloured, clear solution without any haziness

Reaction

Reaction of 3% w/v aqueous solution at 25°C. pH : 7.6±0.2

pH range

7.40-7.80

Cultural Response/Characteristics

DM1071: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.



Dehydrated Culture Media
Bases / Media Supplements

Organism	Inoculum (CFU)	Growth
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	Luxuriant
<i>Salmonella Typhi</i> ATCC 6539	50-100	Luxuriant
<i>Staphylococcus aureus</i> ATCC 25923	$\geq 10^3$	luxuriant

Storage and Shelf Life

Dried media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8⁰ in sealable plastic bags for 2-5 days.

Further Reading

1. Colle, J.G., Duguid J.P., Fraser A.G. and Marmion, B.P. (Eds.) 1989 Mackie and McCartney Practical Medical Microbiology, Vol. 2, p:134 Longman Group, UK.
2. Watson, K.C. 1955, Isolation of Salmonella Typhi from the blood stream. J. of „Lab and Clinical Medicine 46:128-134.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Donot use the products if it fails to meet specificatons for identity and performens parameters.

