

Technical Information

Fungal Agar (Mycological Agar)

Product Code: DM 1094

Application: Fungal Agar (Mycological agar) is recommended for cultivation of fungi. Also used for carrying stock cultures of fungi and for chlamydospores production.

Composition**

Gms / Litre
10.000
10.000
15.000
7.0±0.2

Principle & Interpretation

Earlier media for fungi generally relied on an acidic pH to make the media less suitable for the growth of many bacteria ⁽¹⁾. Fungal Agar is prepared according to the formulation devised by Huppert and Walker ⁽²⁾.

Mycological media are basal media to which antifungal agents may be added for checking their effect on fungi or bacteria to render them selective for isolation and cultivation of fungi. Mycological Agar is used while working with pathogenic fungi.

Papaic digest of soyabean meal in the medium provides nitrogen, vitamins and minerals necessary to support bacterial growth. Dextrose is a carbon source required for the growth of fungi. The pH may be adjusted to 4.0 after autoclaving by adding sterile 10% lactic acid sodium (MS2095)/acetic acid and used for determining yeast and mould counts of carbonated beverages and food products (3).

Methodology

Suspend 35 grams of powder in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates. For preparing selective media, acidify the media upto pH 3.0-4.0 by the addition of two vials of 10% Lactic Acid Solution (MS2095).

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light amber coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 3.5% w/v aqueous solution at 25°C. pH: 7.0±0.2

pH Range 6.80-7.20

Cultural Response/Characteristics

DM 1094: Cultural characteristics observed after an incubation at 25 - 30°C for 48 - 72 hours (For Trichophyton species longer incubation may be required for upto 7days).





Organism	Inoculum (CFU)	Growth	Recovery
Aspergillus brasiliensis ATCC 16404	50-100	luxuriant	
Candida albi cans ATCC 10231	50-100	luxuriant	>=70%
Lactobacillus acidophilus ATCC 11506	50-100	luxuriant	>=70%
Saccharomyces cerevisiae ATCC 9763	50-100	luxuriant	>=70%
Saccharomyces uvarum ATCC 28098	50-100	luxuriant	>=70%
Staphylococcus aureus ATCC 25923	50-100	luxuriant	>=70%
Trichophyton mentagrophytes ATCC 9533	50-100	luxuriant	

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label. **Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

Further Reading

- 1. A. J. Clin. Path., 1951, 21: 684.
- 2. Huppert M., and Walker L. J., 1958, Am. J. Clin. Pathol., 29:291
- 3. Speck M. L., (Eds.), 1984, Compendium of Methods for the Microbiological Examination of Foods, 2nd Ed., APHA, Washington, D.C.

Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performens parameters.

