

## Technical Information

### Micro Vitamin Test Culture Agar

#### Product Code: DM 1132

**Application:** - Micro Vitamin Test Culture Agar is recommended for cultivation and maintenance of stock cultures of *Lactobacilli* used in microbiological assays of vitamins.

#### Composition\*\*

Ingredients	Gms / Litre
Yeast extract	20.000
Peptic digest of animal tissue	5.000
Dextrose	10.000
Monopotassium phosphate	2.000
Polysorbate 80	0.100
Agar	15.000
Final pH (25°C)	6.7±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

*Lactobacillus* is a genus of gram-positive facultative anaerobic lactic acid bacteria. The lactic acid bacteria are so named because most of its members convert lactose and other sugars to lactic acid. They are common and usually benign. Many species are prominent in decaying plant material. The production of lactic acid makes its environment acidic which inhibits the growth of some harmful bacteria. Three types of media are generally used in microbiological assays namely maintenance media, inoculum /cultivation media and the test assay media.

Micro Vitamin Test Agar is used for maintaining stock cultures of and other test organisms used in microbiological assays <sup>(1)</sup>. This media can also be used for routine cultivation of in microbiological assays of vitamins and in inoculum preparation for assays.

Peptic digest of animal tissue and yeast extract in the medium provide nitrogen, sulphur, vitamins and other essential nutrients for growth. Dextrose is the energy source. Polysorbate 80 is the fatty acid source. Monopotassium phosphate buffers the medium.

Stock cultures are prepared by stab inoculation in triplicates. One is used for preparation of stock cultures while others are used for inoculum preparation for assays. Transfer of cultures should be made at weekly or biweekly intervals.

Suspend a 16-24 hours culture of from Micro Vitamin Test Culture Agar into Micro Vitamin Test Inoculum Broth. After an incubation at 35-37°C for 18-24 hours, centrifuge the culture and decant the supernatant. Re-suspend the centrifuged cells in 10 ml of sterile saline suspension. Repeat the washing two more times. Dilute the washed cell suspension with basal assay medium or as desired to obtain the required density of cells.

For procedure of Vitamin Assay, refer standard references <sup>(2)</sup>.

#### Methodology

Suspend 52.1 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Dispense and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

#### Quality Control

##### Physical Appearance

Cream to yellow homogeneous free flowing powder

##### Gelling

Firm, comparable with 1.5% Agar gel

**Colour and Clarity of prepared medium**

Light yellow coloured clear to slightly opalescent gel forms in tubes as butts

**Reaction**

Reaction of 5.2% w/v aqueous solution at 25°C pH : 6.7±0.2

**pH range** 6.50-6.90

**Cultural Response/ characteristics**

DM1132: Cultural characteristics observed after an incubationt 35-37°C for 24 - 48 hours.

Organism	Inoculum (CFU)	Growth
Lactobacillus casei ATCC 9595	50-100	good-luxuriant
Lactobacillus viridescens ATCC 12706	50-100	good-luxuriant
Lactobacillus leichmanni ATCC 4797	50-100	good-luxuriant
Lactobacillus plantarum ATCC 8014	50-100	good-luxuriant

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. Atlas R. M., 1993, Handbook of Microbiological Media, Parks L.C., (Ed.), CRC Press, Inc.
2. Horwitz, (Ed.), 2000, Official Methods of Analysis of AOAC International, 17th Ed., Vol. I, AOAC International, Gaithersburg, Md.

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specificatons for identity and performens parameters.