

## Technical Information

### Antibiotic Assay Medium No.9 (Polymyxin Base Agar)

#### Product Code: DM1147

**Application:** - Antibiotic Assay Medium No.9 (Polymyxin Base Agar) is used as base layer for plate assay of Carbenicillin, Colistin, Colistimethate sodium and Polymyxin B.

#### Composition\*\*

Ingredients	Gms / Litre
Casein enzymic hydrolysate	17.000
Papaic digest of soyabean meal	3.000
Sodium chloride	5.000
Dextrose	2.500
Dipotassium phosphate	2.500
Agar	20.000
Final pH ( at 25°C)	7.2±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

This medium is widely recommended for assay of Polymyxin B, Colistimethate sodium and Colistin using *Bordetella bronchiseptica* as test organisms. Carbenicillin assay is also performed using this medium with *Pseudomonas aeruginosa*. The medium is numerically identical with the name assigned by Groove and Randall <sup>(1)</sup>.

Casein enzymic hydrolysate and papaic digest of soyabean meal serves as source for essential nutrients. Dextrose stimulates the growth by providing carbon and energy. Phosphates in the medium enhance buffering action and sodium chloride maintains osmotic equilibrium in the medium. Agar concentration provides control over the diffusion activity of polymyxin B antibiotics and provides solid substratum to support the seed agar layer.

To perform the antibiotic assay the Base Agar should be prepared on the same day as the test. For the cylinder method, a base layer of 21 ml is required. Once the base medium has solidified, seed layer inoculated with the standardized culture can be overlaid. Even distribution of the layer is important.

#### Methodology

Suspend 50 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

#### Quality Control

##### Physical Appearance

Cream to yellow homogeneous free flowing powder

##### Gelling

Firm, comparable with 2.0% agar gel.

##### Colour and Clarity of prepared medium

Light amber coloured clear to slightly opalescent gel forms in Petri plates

##### Reaction

Reaction of 5% w/v aqueous solution at 25°C. pH : 7.2±0.2



Dehydrated Culture Media  
Bases / Media Supplements

pH range 7.00-7.40

#### Cultural Response/ characteristics

DM 1147: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth	Recovery	antibiotics assayed
<i>Bordetella bronchiseptica</i> ATCC 4617	50-100	good	50-70%	Colistimethate sodium, Colistin, Polymyxin B
<i>Pseudomonas aeruginosa</i> ATCC 25619	50-100	Luxuriant	>=70%	Carbenicillin
<i>Pseudomonas aeruginosa</i> ATCC 27853	50-100	luxuriant	>=70%	Carbenicillin

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. Grove and Randall, 1955, Assay Methods of Antibiotics Medical Encyclopedia, Inc. New York

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.

