

Technical Information

Malt Agar

Product Code: DM 1253

Application: - Malt Agar is used for the detection and isolation of yeasts and moulds from dairy products, foods and other materials. Also used for carrying stock cultures of yeasts and moulds.

Composition**

Ingredients	Gms / Litre
Malt extract	30.000
Agar	15.000
Final pH (25°C)	5.5±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Media based on malt extract may be considered as general purpose media and is very suitable for the cultivation of fastidious microorganisms. With acidic pH, the media is used for the isolation, cultivation and maintenance of yeast and moulds.

Malt media for yeasts and moulds have been widely used for many years. In 1919, Reddish ⁽¹⁾ prepared a suitable substitute for beer wort from malt extract. Malt Agar is included in Official Methods of Analysis of AOAC International ⁽²⁾. It is recommended by APHA ⁽³⁾ for use in both antibiotic and acidified standard methods for yeast and mould counts in food. This medium is also used for maintaining stock cultures of fungi.

Malt Agar contains malt extract, which provides carbon, protein and nutrient sources required for the growth of microorganisms. The acidified medium inhibits the growth of bacteria and allows good recovery of yeasts and moulds ⁽⁴⁾. Heating process during rehydration and sterilization should be for shorter period as excessive heat causes partial hydrolysis of the agar, which results in inability to gel properly when cooled. If desired additional 5 grams of agar may be added.

Methodology

Suspend 45 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at (118°C) for 15 minutes. Avoid overheating, as it will result in a softer and darker agar.

Quality Control

Physical Appearance

Cream to brownish yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Light amber coloured clear to slightly opalescent gel forms in tubes or Petri plates

Reaction

Reaction of 4.5% w/v aqueous solution at 25°C. pH:-5.5±0.2

pH range 5.30-5.70

Cultural Response/ characteristics

DM 1253: Cultural characteristics was observed after an incubation at 25 - 30°C for 40 - 48 hours.

Organism	Inoculum (CFU)	Growth	Recovery
*Aspergillus brasiliensis ATCC 16404	50-100	luxuriant	
Candida albicans ATCC 10231	50-100	luxuriant	>=70%
Saccharomyces cerevisiae ATCC 9763	50-100	luxuriant	>=70%



Dehydrated Culture Media
Bases / Media Supplements

*Key: Formerly known as *Aspergillus niger* ATCC 16404

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Reddish, 1919, Abstr. Bacteriol., 3:6.
2. Williams, (Ed.), 2005, Official Methods of Analysis of the Association of Official Analytical Chemists, 19th Ed., AOAC, Washington, D.C.
3. Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.
4. Can. Dept. Agr. Pamphlet, 92-N.S.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.

