

Dehydrated Culture Media Bases / Media Supplements

Technical Information

Emerson Agar

Product Code: DM1325

Application: - Emerson Agar is used for isolation and cultivation of Actinomycetaceae, Streptomycetaceae, fungi and moulds.

| Composition** | | |
|--|-----------------|--|
| Ingredients | Gms / Litre | |
| Beef extract | 4.000 | |
| Yeast extract | 1.000 | |
| Peptic digest of animal tissue | 4.000 | |
| Dextrose | 10.000 | |
| Sodium chloride | 2.500 | |
| Agar | 20.000 | |
| Final pH (at 25°C) | 7.0±0.2 | |
| **Formula adjusted standardized to suit perform: | ance parameters | |

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Emerson Agar which was originally formulated by Emerson et ^{al (1)} and used for the cultivation of moulds and bacterial species resembling moulds ⁽³⁾ was further modified by Gottlieb et al ⁽²⁾ for screening potent antibiotic-producing organisms ^{(4).} In their study, they stored Streptomyces in soil for long time and transferred them to slants of Emerson Agar as and when needed. The slant cultures were incubated for 3-7 days. The spores were gently scraped from the cultures surface to form a spore inoculum.

Yeast extract provides a source of trace elements, vitamins and amino acids. For the selective isolation of Streptomyces species, cycloheximide is

incorporated in the medium, which limits the growth of moulds. This medium is also used for routine cultivation and maintenance of pure cultures.

Methodology

Suspend 41.5 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Add 0.05 grams / litre cycloheximide, if desired. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.

Quality Control

Physical Applearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 2.0% agar gel.

Colour and Clarity of prepared medium

Light amber coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 4.15% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH range 6.80-7.20

Cultural Response/ characteristic

DM 1325 Cultural characteristic observed when sub cultured on Emerson Agar at 22-25⁰ C for 2-5 days.





Dehydrated Culture Media Bases / Media Supplements

| Organism |
|--|
| *Aspergillus brasiliensis ATCC 16404 |
| Saccharomyces cerevisiae ATCC 9763 |
| Streptomyces albus subspalbus ATCC 3004 |
| Streptomyces lavendulae ATCC 8664 |
| Streptomyces achromogenes ATCC 12767 |
| *Key: Formerly known as Aspergillus niger ATCC 16404 |

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label. **Prepared Media**: 2-8° in sealable plastic bags for 2-5 days.

Growth luxuriant luxuriant luxuriant luxuriant luxuriant

Further Reading

- 1. Emerson R. L., Whiffen A. J., Bohonos M. and DeBoer C., 1946, J. Bacteriol., 52:357.
- 2. Gottlieb D., Bhattacharya P. K., Anderson H. W. and Carter H. E., 1948, J. Bacteriol., 55:409.
- 3. Haynes W. C., Wickerham L. J. and Hesseltine C. W., 1955, Appl. Microbiol., 3:361.
- 4. Schmitz H. and Woodside R., 1955, Antibiot. Chemother., 5:652.

Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specificatons for identity and performens parameters.

