

Technical Information

Differential Reinforced Clostridial Broth Base

Product Code: DM 1549

Application: Differential Reinforced Clostridial Medium is used for the cultivation of Clostridia from water.

Composition**

Composition		
Ingredients	Gms / Litre	
Peptic digest of animal tissue	10.000	
Beef extract	10.000	
Yeast extract	1.500	
Starch	1.000	
Sodium acetate, hydrated	5.000	
Glucose	1.000	
L-Cysteine hydrochloride	0.500	
Final pH (at 25°C)	7.2±0.2	

Principle & Interpretation

Hirsch and Grinstead ⁽¹⁾ devised Differential Reinforced Clostridial Medium to initiate the growth from small inoculum and get a higher Clostridial count. Later, Barnes and Ingram ⁽²⁾ used the medium to develop vegetative cells in assays of *Clostridium perfringens*. This medium is developed for the isolation of sulphite-reducing Clostridia from food and for their enumeration in water by multiple tube method. Differential Reinforced Clostridial Broth is used to determine the count of sulphite reducing bacteria by MPN technique ⁽³⁾.

Peptic digest of animal tissue, beef extract, yeast extract, starch, sodium acetate provide essential nutrients for bacterial metabolism. Glucose is the fermentable carbohydrate and serves as carbon and energy source. L-cysteine hydrochloride acts as reducing agent. Sodium sulphite and ferric citrate are added as indicators. Sulphite reducing clostridia produce sulphide from sulphite, which results in the formation of black coloured medium.

Methodology

Suspend 29 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Just before use add 0.5 ml filter sterilized solution, prepared by mixing equal volumes of 4% w/v solution of sodium sulphite and 7% w/v ferric citrate, to 25 ml of single strength medium or 0.4 ml and 2 ml to 10 ml and 50 ml of double strength medium respectively. Mix well.





Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Light yellow coloured, clear solution without any precipitate

Reaction

Reaction of 2.9% w/v aqueous solution at 25°C. pH: 7.2±0.2

pH range

7.00-7.40

Cultural Response/Characteristics

DM 1549: Cultural characteristics observed in an anaerobic atmosphere, with added 4% w/v solution of Sodium sulphite and 7% w/v Ferric citrate after an incubation at 30-35°C within 1 week.

Organism	Inoculum (CFU)	Growth	H2S production
Clostridium perfringens ATCC 13124	50-100	Good-luxuriant	positive reaction, blackening of medium
Clostridium sporogenes ATCC 11437	50-100	Good-luxuriant	positive reaction, blackening of medium

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

- 1. Hirsch A. and Grinstead E., 1954, J. Dairy Res. 21:101
- 2. Barnes E. M. and Ingram M., 1956, J. Appl. Bacteriol., 19(1):117.
- 3. Collee J. G., Fraser A. G., Marmion B. P., Simmons A., (Eds.), Mackie and McCartney, Practical Medical Microbiology, 1996, 14th Edition, Churchill Livingstone

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performens parameters.

