

Technical Information

MY 40 Agar (Osmophilic Agar)

Product Code: DM 1594

Application: - MY 40 Agar (Osmophilic Agar) is used for detection and isolation of osmophilic microorganisms from food samples.

Composition**

Ingredients	Gms / Litre
Malt extract	20.000
Yeast extract	5.000
Sucrose	400.000
Agar	20.000

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Osmophilic yeasts usually are responsible for the spoilage of food with high sugar content including jams, honey, concentrated fruit juices, chocolate candy with soft centres etc. ^(1, 2). Though osmophilic yeasts are of no public health significance, yet they are of economic importance to the food industry ⁽²⁾. A simple presence-absence test for detection of small numbers of osmotolerant yeasts in high-sugar foods is useful for enumeration. Osmophilic Agar (MY 40 Agar) is generally used for this purpose. Walker and Ayers ⁽³⁾ in their review have differentiated between osmophilic yeasts and osmoduric yeasts. Almost all of the known osmophilic yeasts are *Saccharomyces* species e.g. *Saccharomyces rouxii, var. polymorphus, Saccharomyces mellis* etc. Improved recovery of osmophillic yeasts has been reported on media, which resemble the composition of the food under examination or contain high sugar concentrations ⁽⁴⁾.

Osmophilic Agar is recommended for cultivation of a wide variety of osmophilic organisms ⁽⁵⁾. MY in MY-40 Agar stands for malt extract and yeast extract and 40 for the 40% of sucrose in the medium, which meets the above requirements.

The medium contains malt extract and yeast extract which supply the nitrogenous nutrients, amino acids, vitamins, trace ingredients to the osmophilic yeasts. 40% sucrose in the medium satisfy the nutritional need of these yeasts.

Methodology

Suspend 44.5 grams of powder media in 100 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. DO NOT OVERHEAT. Mix well and pour into sterile Petri plates.

Quality Control

Physical Appearance

Off-white to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 2.0% agar gel.

Colour and Clarity of prepared medium

Medium amber coloured slightly opalescent gel forms in Petri plates

Cultural Response/ characteristices

DM 1594: Cultural characteristics observed after an incubation at 25-30°C for upto 1 week.





Organism	Growth
*Aspergillus brasiliensis ATCC 16404	good-luxuriant
Mucor racemosus ATCC 22365	good-luxuriant
Pencillium notatum ATCC 10108	good-luxuriant
Saccharomyces rouxii ATCC 28253	good-luxuriant

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

- 1. Rose A. H. and Harrison J. S., (Eds.), 1970, The Yeasts, Vol. 3, Academic Press, New York.
- 2. Tilbury R. H., 1980, "Biology and Activities of Yeasts", Skinner and others (Ed.), Academic Press, London.
- 3. Walker H.W., and Ayers J.C., 1970, Jn. A.H. Rose and J.S. Harrison (Eds.), `The Yeasts', Vol. 3, Academic Press, Inc., New York.
- 4. Anand J.C. and Brown A.D., 1968, J. Gen. Microbiol., 52:205.
- 5. Atlas R. M., 2004, Handbook of Microbiological Media, Lawrence C. Parks (Ed.), 3rd Edition, CRC Press.

Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performens parameters.

