



Dehydrated Culture Media  
Bases / Media Supplements

## Technical Information

### Esculin Azide Broth

**Product Code: DM 1749**

**Application:** - Esculin Azide Broth is used for selective cultivation and identification of Streptococci.

### Composition\*\*

Ingredients	Gms / Litre
Peptic digest of animal tissue	20.000
Yeast extract	5.000
Bile salts	10.000
Sodium citrate	1.000
Esculin	1.000
Ferric ammonium citrate	0.500
Sodium azide	0.250
Final pH ( at 25°C)	7.2±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Principle & Interpretation

Esculin Azide Broth is prepared based on the modification of original formula of Isenberg in which bile concentration was 40 g/l. In Esculin Azide Broth, the concentration of bile was reduced to 10 g/l and also additional sodium azide was added <sup>(1)</sup>. Enterococci are able to hydrolyze esculin whereas other streptococci are not able to do so. Swan <sup>(2)</sup> used an esculin medium containing 40% bile salts and reported that a positive reaction on bile esculin medium could be correlated with a serological group D precipitin reaction. Further studies by Facklam and Moody presumptively identified group D Streptococci and found that the bile esculin test provided a reliable means of identifying group D Streptococci and differentiating them from non-group D Streptococci <sup>(3)</sup>

The broth is selective due to presence of bile salts and sodium azide and provides rapid growth of Streptococci. Peptic digest of animal tissue and yeast extract provide nitrogenous nutrients to the organisms. Bile salts inhibit other gram-positive bacteria while sodium azide inhibits gram-negative bacteria. Streptococci hydrolyze esculin to esculetin and dextrose. Esculetin and ferric ammonium citrate forms dark brown to black complex, imparting dark brown colour to the broth.

### Methodology

Suspend 37.75 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Dispense as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Warning: Sodium azide has a tendency to form explosive metal azides with plumbing materials. It is advisable to use enough water to flush off the disposables.

### Quality Control

#### Physical Appearance

Cream to yellow homogeneous free flowing powder

#### Colour and Clarity of prepared medium

Amber coloured, clear solution having slight purplish tinge

#### Reaction

Reaction of 3.78% w/v aqueous solution at 25°C. pH : 7.2±0.2

pH range 7.00-7.40



### Cultural Response/ characteristics

DM 1749: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth	Esculin Hydrolysis
<i>Enterococcus faecalis</i> ATCC 29212	50-100	good to luxuriant	Positive reaction, blackening of medium
<i>Escherichia coli</i> ATCC 25922	$\geq 10^3$	inhibited	
<i>Streptococcus bovis</i> ATCC 27960	50-100	good to luxuriant	Positive reaction, blackening of medium
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	poor	Negative reaction, no change
<i>Proteus mirabilis</i> ATCC 25933	50-100	poor to fair	Negative reaction, no change

### Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

### Further Reading

1. Isenberg, 1970, Clin. Lab. Forum. 2. Swan A., 1954, J. Clin. Pathol., 7:1603. Facklam R. R. and Moody M. D., 1970, Appl. Microbiol., 20:2451. Isenberg, 1970, Clin. Lab. Forum.
2. Swan A., 1954, J. Clin. Pathol., 7:160
3. Facklam R. R. and Moody M. D., 1970, Appl. Microbiol., 20:245

### Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performs parameters.