



Dehydrated Culture Media
Bases / Media Supplements

STechnical Information

Skim Milk Agar

Product Code: DM 1763

Application: - Skim Milk Agar is used for cultivation and enumeration of microorganisms encountered in dairy industry.

Composition**

Ingredients	Gms / Litre
Skim milk powder	28.000
Casein enzymic hydrolysate	5.000
Yeast extract	2.500
Dextrose	1.000
Agar	15.000
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Skim Milk Agar is used for the demonstration of coagulation and proteolysis of casein ⁽¹⁾. The medium is recommended by APHA ⁽²⁾ for cultivation and enumeration of microorganisms encountered in dairy industry ⁽³⁾. Addition of skim milk to any nutrient-rich medium creates favourable conditions for growth of organisms, which are present in milk. The number of bacteria isolated on this media thus is more than the number of organisms isolated on a regular medium ⁽⁴⁾. Proteolytic bacteria hydrolyze casein to form soluble nitrogenous compounds indicated as clear zone surrounding the colonies. More clear zones are seen on milk agar if, the bacteria produce acid from fermentable carbohydrates in the medium. Casein enzymic hydrolysate provides amino acids and other complex nitrogenous substances. Yeast extracts supplies vitamin B complex. Addition of skim milk in the medium makes the conditions optimal for microorganisms encountered in milk. Glucose acts as the carbon source.

Methodology

Suspend 51.5 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Off white coloured opaque gel forms in Petri plates

Reaction

Reaction of 5.15% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH Range:-

6.80-7.20

Cultural Response/Characteristics

DM 1763: Cultural characteristics observed after anincubation at 35-37°C for 18-24 hours.





Dehydrated Culture Media
Bases / Media Supplements

Organism	Inoculum(CFU)	Growth	Recovery	Colour of Colony
<i>Bacillus subtilis</i> ATCC 6633	50-100	luxuriant	>=70%	negative reaction, no clear zone surrounding colonies
<i>Enterococcus faecalis</i> ATCC 29212	50-100	Good- luxuriant	>=70%	negative reaction, no clear zone surrounding colonies
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	>=70%	positive reaction, clear zone surrounding colonies
<i>Proteus mirabilis</i> ATCC 25933	50-100	luxuriant	>=70%	positive reaction, clear zone surrounding colonies
<i>Pseudomonas aeruginosa</i> ATCC 27853	50-100	luxuriant	>=70%	positive reaction, clear zone surrounding colonies
<i>Serratia marcescens</i> ATCC 8100	50-100	luxuriant	>=70%	positive reaction, clear zone surrounding colonies

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Frazier W. C. and Ripp P., 1928, J. Bacteriol., 16: 57.
2. Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C.
3. Wehr H. M. and Frank J. H., 2004, Standard Methods for the Microbiological Examination of Dairy Products, 17th Ed., APHA Inc., Washington, D.C.
4. Terplan G. Rundfeldt, H.u. Zaadhof, K.J. Zur Eignung verschiedener Nährböden für die Bestimmung der Gesamtkeimzahl der Milch. - Arch. Lebensmittelhyg., 18; 9-11 (1967).

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performance parameters.

