

# **Technical Information**

### **Antibiotic Assay Medium M- AOAC**

Product Code: DM 1992

**Application:** - Antibiotic Assay Medium M is recommended for microbiological assay of Lasalocid using *Bacillus subtilis* ATCC 6633 as test organism.

Composition\*\*

Composition		
Ingredients	Gms / Litre	
Yeast extract	2.500	
Dextrose	10.000	
Dipotassium hydrogen phosphate	0.690	
Potassium dihydrogen phosphate	0.450	
Agar	20.000	
Final pH ( at 25°C)	$6.0 \pm 0.2$	
**Formula adjusted, standardized to suit perforn	nance	
parameters		

**Principle & Interpretation** 

Antibiotic Assay Medium M is formulated in accordance with AOAC (1) for the microbiological assay of Lasalocid in feeds, using Bacillus subtillis (ATCC 6633) as the test organism.

Prepare slant culture of *Bacillus subtillis* (ATCC 6633) on Assay Medium No. 1 and incubate for 16-24 hours at 37°C. Wash the growth with sterile distilled water and transfer it to surface of Assay Medium No. 32 and incubate at 37°C for 7 days. Wash the growth with sterile distilled water. Heat to 65°C for 30 minutes in water bath. Centrifuge, decant the supernatant and resuspend the cells. Repeat this for 3 minutes in water bath. Dilute suspension with sterile distilled water (1 + 50) to read 20%T on sprectrophotometer at 530 nm before use.

Use single inoculated agar layer. For this optimum concentration of suspension of *Bacillus subtillis* is determined prior to assay to be added to Medium M to obtain inhibition zone of adequate size (17.5 ± 2.5 mm with 1.0 µg/ml). For actual assay add appropriate amount of suspension to sterile, molten medium M (pH 6.0). Mix and add 6 ml to each plate. Prepare plates two to three hours before use. Weigh 1.0 g premix. Transfer to flask and add 100 ml methanol. Shake vigorously for 3 minutes and dilute 4 ml of this to 100 ml methanol. Further dilute 3 ml with 22 ml methanol and water to 100 ml (1 ml= ca/µg lasalocid Na/ml 25% methanol). Prepare final concentration of feed to 0.0075%. For more details refer AOAC.

Using lasalocid, sodium obtains standard response line, assay solution. Place cylinders on each plate and alternatively fill with reference concentration and other standard concentration. Incubate at 35-36°C. Calculate zone diameters of L (Low concentration giving measurable zone) and H (Highest concentration) of standard response line and connect with straight line. This corrected reference point is used for sample calculations. Average the 9 readings of reference concentration and 9 readings of assay solution. If assay solution gives larger average than reference concentration, add difference between them to reference point on standard response line. If the assay solution gives lower average than reference concentration, subtract the difference from reference point. Using the corrected value of assay solution, amount of antibiotic is determined.

### Methodology

Suspend 33.64 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well before pouring in sterile Petri plates.





## **Quality Control**

### **Physical Appearance**

Cream to yellow homogeneous free flowing powder

#### Gelling

Firm, comparable with 2.0% Agar gel

### Colour and Clarity

Yellow coloured clear to slightly opalescent gel forms in Petri plates

#### Reaction

Reaction of 3.3 6% w/v aqueous solution at 25°C. pH: 6.0±0.2

pH range 5.80-6.20

#### Cultural Response/ characteristices

DM 1992: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.

Organism Growth zones with Cultural Response Inhibition

Bacillus subtilis ATCC 6633 luxuriant Lasalocid

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label. **Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

### **Further Reading**

1. Williams (Ed.), 2005, Official Methods of Analysis of AOAC International, 19th ed., AOAC, International, Washington D. C.

### Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performens parameters.

