



Dehydrated Culture Media  
Bases / Media Supplements

## Technical Information

### Esculin Iron Agar

#### Product Code: DM 2044

**Application:** - Esculin Iron Agar is used for verification of enterococcal colonies on membrane filters through which water samples have been filtered and which have been incubated on M-Enterococcus Agar, Modified (M1048).

#### Composition\*\*

Ingredients	Gms / Litre
Esculin	1.000
Ferric ammonium citrate	0.500
Agar	15.000
Final pH ( at 25°C)	7.1±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

Enterococci are indicators for the bacteriological quality of recreational waters, since they occur in faeces of humans and warm-blooded animals <sup>(1)</sup>. Detection and quantitation of Enterococci is necessary because gastroenteritis is associated with swimming in recreational water, which in turn is related to enterococcal densities <sup>(2)</sup>. Esculin Iron Agar is used in conjunction with M - Enterococcus Agar, Modified, (DM2048) for verification of enterococcal colonies in fresh and marine recreational water, which is also followed by APHA <sup>(3)</sup>. Esculin in the medium is hydrolyzed by Enterococci to form esculetin and dextrose. Esculetin reacts with the iron salt (ferric ammonium citrate) and produces a dark brown to black complex, which appears around the colonies.

In the membrane filtration technique, two media, namely M-Enterococcus Agar, Modified (DM2048) and Esculin Iron Agar (DM2044) are used in conjunction, where the former serves as a selective medium while the later confirms the identification of colonies on the basis of its ability to hydrolyze esculin. The membrane filter used to filter the test water sample is aseptically placed on M-Enterococcus Agar, Modified (DM2048) and incubated at 40-42°C for 48 hours. After incubation the membrane is aseptically transferred to Esculin Iron Agar (DM2044) plate and incubated at 40-42°C for 20 minutes. After incubation count and record the number of pink to red colonies with black or reddish brown precipitate on the underside of the membrane. If required, magnifying glass or fluorescent lamp may be used for counting the visible colonies. Following formula is used for the final calculation <sup>(3)</sup>.

Enterococci / 100 ml = No of enterococcal colonies X 100

Volume of sample filtered

#### Methodology

Suspend 16.5 grams of powder media in 1000 ml distilled water. Shake well & heat to boiling with frequent stirring. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45°C and pour into sterile Petri plates to a depth of 4-5 mm.

#### Quality Control

##### Physical Appearance

Light yellow to light brown homogeneous free flowing powder

##### Gelling

Firm, comparable with 1.5% Agar gel

##### Colour and Clarity of prepared medium

Medium amber coloured, clear to slightly opalescent gel forms in Petri plates

##### Reaction

Reaction of 1.65% w/v aqueous solution at 25°C. pH : 7.1±0.2

pH range 6.90-7.30





Dehydrated Culture Media  
Bases / Media Supplements

#### Cultural Response/ characteristics

DM 2044: Cultural characteristics observed after an incubation at 40-42°C for 18-24 hours on M-Enterococcus Agar, Modified (DM1048) and after 20 minutes at 40-42°C on Esculin Iron Agar (MS2044).

Organism	Growth	Esculin Hydrolysis	Inoculum(CFU)	Colour of colony
Escherichia coli ATCC 25922	none-poor	negative reaction		
Enterococcus faecalis ATCC 29212	good-luxuriant	positive reaction, brown to black precipitate round colonies.	50-100	pink to red

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. U. S. Environmental Protection Agency, 1997, EPA Method 1600: Membrane Filter Test Method for Enterococci in Water, EPA-821-R-97-004, Washington, D.C.
2. Cabelli et al, 1979, Am. J. Public Health, 69:690.
3. Eaton A. D., Clesceri L. S. and Greenberg A. E., (Ed.), 1998, Standard Methods for the Examination of Water and Wastewater, 20th Ed., American Public Health Association, Washington, D.C.

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.

