



Dehydrated Culture Media  
Bases / Media Supplements

## Technical Information

### Tryptone Tellurite Agar Base

**Product Code: DM 2056**

**Application:** - Tryptone Tellurite Agar Base with supplements is used for selective isolation of pathogenic microorganisms from clinical specimens especially from nose, throat and vagina.

### Composition\*\*

Ingredients	Gms / Litre
Casein enzymic hydrolysate	10.000
Peptic digest of animal tissue	10.000
Sodium chloride	5.000
Dextrose	2.000
Agar	20.000
Final pH ( at 25°C)	7.5±0.2

\*\*Formula adjusted, standardized to suit performance parameters

### Principle & Interpretation

Tryptone Tellurite Agar was prepared according to the specifications of the Maryland State, Department of Health as a medium for isolation of members of the genus *Corynebacterium*, particularly in the laboratory diagnosis of diphtheria<sup>(1)</sup>. Pathogenic organisms like *Corynebacterium*, Streptococci are often responsible for infections in throat and nasal cavity<sup>(2)</sup>. Isolation of these organisms is important for diagnosis of fatal infections so that timely action taken may save the life of patient.

Casein enzymic hydrolysate and peptic digest of animal tissue present in the medium provide nitrogenous, carbonaceous compound and trace elements. Dextrose is an energy source. Sodium chloride maintains osmotic equilibrium. Potassium tellurite inhibits gram-negative and most gram-positive bacteria except *Corynebacteria* and some Streptococci<sup>(3)</sup>. Serum provides additional growth factors for *Corynebacterium* species.

### Methodology

Suspend 47 grams of powder media in 950 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 118-121°C for 15 minutes. Cool to around 50-55°C and aseptically add sterile, 10 ml of 1% Potassium Tellurite Solution (MS2052) and 50 ml of sterile sheep, rabbit or human serum. Mix well and pour into sterile Petri plates.

### Quality Control

#### Physical Appearance

Cream to yellow homogeneous free flowing powder

#### Gelling

Firm, comparable with 2.0% Agar gel.

#### Colour and Clarity of prepared medium

Yellow coloured opalescent gel forms in Petri plates.

#### Reaction

Reaction of 4.7% aqueous solution at 25°C. pH : 7.5±0.2

#### pH Range:-

7.30-7.70

#### Cultural Response/Characteristics

DM 2056: Cultural characteristics observed after an incubation at 35- 37°C with added Potassium Tellurite Solution (FD052) and sterile human serum for 24-48 hours.





Dehydrated Culture Media  
Bases / Media Supplements

Organism	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Corynebacterium diphtheria</i> ATCC 11913	50-100	luxuriant	>=50%	grey-black with dark centers
<i>Escherichia coli</i> ATCC 25922	>=10 <sup>8</sup>	inhibited	0%	
<i>Staphylococcus aureus</i> ATCC 25923	50-100	none-poor	<=10%	white-grey

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8<sup>0</sup> in sealable plastic bags for 2-5 days.

## Further Reading

1. Albers W. D., 1947, U.S. Naval Med. Bull., 47:33.
2. Collee J. G., Fraser A. G., Marmion B. P., Simmons A. (Eds.), Mackie and McCartney, Practical Medical Microbiology, 1996, 14th Edition, Churchill Livingstone.
3. MacFaddin J. F., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performance parameters.

