

Technical Information

M17 Agar with Glycerophosphate

Product Code: DM 2063

Application: - M17 Agar with Glycerophosphate is used for cultivation of lactic Streptococci and plaque assay of lactic bacteriophages.

Composition**

Ingredients	Gms / Litre
Papaic digest of soyaben meal	5.000
Biopeptone	5.000
Yeast extract	2.500
Beef extract	5.000
Lactose	5.000
Ascorbic acid	0.500
Magnesium sulphate	0.250
Disodium - β - glycerophosphate	19.000
Agar	10.000
Final pH (at 25°C)	7.1 \pm 0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

M17 media is devised by Terzaghi and Sandine ⁽¹⁾ for the cultivation and enumeration of lactic Streptococci and their bacteriophages. It is possible to study plaque morphology and lysogeny. M17 Agar is recommended by the International Dairy Federation ⁽⁵⁾ for selective enumeration of *Streptococcus thermophilus* from yoghurt. M17 Agar is recommended by APHA for the cultivation of lactic Streptococci ⁽⁶⁾. Lactic Streptococci are nutritionally fastidious and require complex media for optimal growth ^(2, 3). Disodium glycerophosphate maintains the pH above 5.7. The maintenance of pH is very important as lower pH results in injury and less recovery of lactic Streptococci. Glycerophosphate does not form precipitate with calcium which is needed for the plaque assay of lactic bacteriophages.

Papaic digest of soyabean meal, yeast extract, beef extract and biopeptone provide carbonaceous, nitrogenous compounds, vitamin B complex and other essential growth factors. Lactose is the fermentable carbohydrate and ascorbic acid is stimulatory for the growth of lactic Streptococci. Magnesium sulphate provides essential ions to the organisms. Disodium glycerophosphate suppresses *Lactobacillus bulgaricus*. Shankar and Davies ⁽⁴⁾ reported isolation and enumeration of *Streptococcus thermophilus* from yoghurt. DM1017 Agar is also recommended by the International Dairy Federation ⁽⁵⁾ for selective enumeration of *Streptococcus thermophilus* from yoghurt. It is also suitable for cultivation and maintenance of starter cultures for cheese and yoghurt manufacturing. This medium helps in detecting streptococcal mutants that are lactose non-fermenters.

Suggested technique to enumerate streptococci is to seed in mass or by stabbing with agar, melted and cooled to 50-55°C, and incubating them at 42°C for 24 hours period. With these conditions, all the colonies might be streptococci. Longer incubation periods or lower temperatures may cause morphological changes in the colonies, which hinders in the recognition of the colonies. Lactose-positive colonies of streptococci are visible after 15 hours and after 5 days they may reach a diameter of about 3-4 mm, whereas those of lactose-negative are 1 mm in diameter. Bacteriophages presence is observed by appearance of characteristic plaques over the bacterial growth.

Methodology

Suspend 52.25 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and dispense as desired.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.0% Agar gel.

Colour and Clarity of prepared medium

Light yellow coloured slightly opalescent gel forms in Petri plates

Reaction

Reaction of 5.23% w/v aqueous solution at 25°C. pH : 7.1±0.2

pH range: 6.90-7.30

Cultural Response/Characteristics

DM2063: Cultural characteristics observed after an incubation at 35-37°C for 24-48 hours.

Organism	Inoculum (CFU)	Growth	Recovery
<i>Enterococcus faecalis</i> ATCC 29212	50-100	Good-luxuriant	>=50%
<i>Lactobacillus bulgaricus</i> ATCC 11842	50-100	None-poor	>=10%
<i>Lactobacillus leichmannii</i> ATCC 4797	50-100	Good-luxuriant	>=50%
<i>Lactobacillus plantarum</i> ATCC 8014	50-100	Good-luxuriant	>=50%
<i>Streptococcus thermophilus</i> ATCC 14485	50-100	Good-luxuriant	>=50%

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Terzaghi B.E. and Sandine W.E., 1975, Appl. Microbiol., 29:807.
2. Anderson A.W. and Elliker P.R., 1953, J. Dairy Sci., 36:161.
3. Reiter B. and Oran J.D., 1962, J. Dairy Res., 29:63.
4. Shankar P.A. and Davies F.L., 1977, Soc. Dairy Technol., 30:28.
5. International Dairy Federation, 1981, Joint IDF/ISO/AOAC Group E44.
6. Downes F. P. and Ito K. (Eds.), 2001, Compendium of Methods for Microbiological Examination of Food, 4th Ed., APHA, Washington, D.C.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performance parameters.