



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

M-(HPC) Heterotrophic Plate Count Agar Base

Product Code: DM 2123

Application: - M-(HPC) Heterotrophic Plate Count Agar Base is recommended for the enumeration of heterotrophic microorganisms from water samples using membrane filter technique.

Composition**

Ingredients	Gms / Litre
Peptic digest of animal tissue	20.000
Gelatin	25.000
Agar	15.000
Final pH (25°C)	7.1±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Heterotrophs are organisms comprising of bacteria, yeasts and moulds that require an external source of organic carbon for growth. The heterotrophic plate count (HPC), formerly known as the standard plate count is a procedure for counting the number of live heterotrophic bacteria in water and measuring changes during water treatment and distribution or in swimming pools. Heterotrophic Plate Count Method has been used in many different forms and is widely used to measure the heterotrophic microorganism population in drinking water and other media. M-(HPC) Heterotrophic Plate Count Agar Base with added glycerol is recommended for the detection of heterotrophic organisms in potable water, swimming pool and other water sources^(1, 2). Three different methods are described for determining the heterotrophic plate count i.e. pour plate method, spread plate method and membrane filter method. M-(HPC) Heterotrophic Plate Count Agar Base is used in the membrane filtration technique. M-(HPC) Heterotrophic Plate Count Broth Base can also be employed for the determination of Heterotrophic Plate Count by the membrane filter method. Sterile cotton absorbent pads are saturated with the broth medium. Membrane filters are then placed on these saturated cotton absorbent pads or agar surface and incubated.

Peptic digest of animal tissue is the source of nutrients for organisms, which are not highly fastidious. Gelatin is utilized by microorganisms through a proteolytic mechanism. The addition of glycerol to the basal medium provides a source of carbon and energy.

Methodology

Suspend 60.0 grams of powder media in 1000 ml distilled water containing 10 ml glycerol. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 5 minutes. Mix well and dispense as desired.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel and 2.5% Gelatin gel

Colour and Clarity of prepared medium

Light yellow coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 6.0% w/v aqueous solution at 25°C. pH : 7.1±0.2

pH range 6.90-7.30

Cultural Response/ characteristics

DM 2123: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.





Dehydrated Culture Media
Bases / Media Supplements

Organism	Inoculum (CFU)	Growth
Escherichia coli ATCC 25922	50-100	luxuriant
Enterococcus faecalis ATCC 29212	50-100	luxuriant
Pseudomonas aeruginosa ATCC 27853	50-100	luxuriant

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Taylor R. H. and Geldreich E. E., 1979, J. Am. Water works Assoc. 71:402.
2. Eaton A. D., Clesceri L. S. and Greenberg A W., (Eds.), 2005, Standard Methods for the Examination of Water and Wastewater, 21st Ed., APHA, Washington, D.C.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Donot use the products if it fails to meet specificatons for identity and performens parameters.

