

Technical Information

Coliform Broth

Product Code: DM 2211

Application: - Coliform Broth is recommended for isolation and cultivation of coliform organisms from cream, yogurt and raw milk

Composition**

Ingredients	Gms / Litre	
Proteose peptone	10.000	
Yeast extract	6.000	
Bile salts	20.000	
Sodium deoxycholate	0.100	
Lactose	20.000	
Sodium lauryl sulphate	1.000	
Bromocresol purple	0.035	
Final pH (at 25°C)	7.0±0.2	
**=		

^{*}Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Bacteriological Examination of water samples to determine is the presence of total coliforms, faecal coliforms or *Escherichia coli* by the most probable number (MPN) procedures or the membrane filter (MF) technique ⁽²⁾ is well recognized indication of unsafe or poor water quality for which corrective measures should be taken. Coliform Broth is recommended for isolation and cultivation of coliforms organisms from cream yogurt and raw milk ⁽¹⁾.

Proteose peptone and yeast extract provides nitrogenous and carbonaceous compounds, vitamin B complex. Lactose is the fermentable carbohydrate. Bromo cresol purple is the pH indicator. Coliforms that ferment lactose produce acid and gas. The acidity formed is indicated by a colour change of the medium from purple to yellow, by the pH indicator dye bromo cresol purple. Sodium deoxycholate and bile salts inhibit gram positive bacteria. Sodium lauryl sulphate is inhibitory to many organisms but not to coliforms. A distinct yellow colour results from the fermentation of lactose and gas production can be detected as bubbles with gentle shaking.

Methodology

Suspend 57.14 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Dispense as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Purple coloured, clear solution without any precipitate

Reaction

Reaction of 5.7 1% w/v aqueous solution at 25°C. pH : 7.0±0.2

pH range 6.80-7.20





Cultural Response/ characteristics

DM 2211: Cultural characteristics observed after an incubation at 35-37°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth	Acid	Gas
Enterobacter aerogenes ATCC 13048	50-100	good-luxuriant	positive reaction, yellow colour	positive reaction
Escherichia coli ATCC 25922	50-100	good-luxuriant	positive reaction, yellow colour	positive reaction
Staphylococcus aureus ATCC 25923	>=10 ³	inhibited		
Salmonella Typhimurium ATCC 14028	50-100	good-luxuriant	negative reaction, no colour change	negative reaction

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

- 1. Atlas R. M., 2004, Handbook of Microbiological Media, Lawrence C. Parks (Ed.), 3rd Edition, CRC Press.
- 2. Eaton A. D., Clesceri L. S., Rice E. W. and Greenberg A. E., (Eds.), 2005, Standard Methods for the Examination of Water and Wastewater, 21st Ed., APHA, Washington, D.C.

Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performens parameters.

