

# **Technical Information**

## Tartoff - Hobbs Broth (Terrific Broth)

## Product Code: DM 2250

Application: - Tartoff-Hobbs (Terrific) Broth is recommended for the cultivation of recombinant strains of Escherichia coli.

Composition**		
Ingredients	Gms / Litre	
Casein enzymic hydrolysate	12.000	
Yeast extract Monopotassium phosphate	24.000 2.200	
Dipotassium phosphate	9.400	
Final pH ( at 25°C) **Formula adjusted, standardized to suit performance p	7.2±0.2 arameters	

## Principle & Interpretation

Terrific Broth Medium was developed by Tartoff Hobbs for cultivation of recombinant *Escherichia coli* strains. These strains have extended growth phase when cultivated in this medium <sup>(1).</sup> Tartoff-Hobbs Medium supports high cellular density and mass and maintains the growth in the logarithmic phase for a long time. Due to this fact, it provides greater yields of recombinant proteins and plasmid DNA. Often, Tartoff-Hobbs (Terrific) Broth substitutes Luria Bertani Broth (DM2151), to get enhanced yields of plasmid DNA and recombinant proteins. The procedures for inoculation, incubation and generation of recombinant strains are detailed by Sambrook et al <sup>(2).</sup>

Casein enzymic hydrolysate and yeast extract supply the necessary nutrients and cofactors for the excellent growth of recombinant strains of *E. coli*. The addition of extra peptone and yeast extract in the medium allows higher plasmid yield per volume. The two phosphates provide good buffering action to the medium. Glycerol is used as the carbohydrate source. Unlike glucose, glycerol is not fermented to acetic acid.

#### Methodology

Suspend 47.6 grams of powder media in 1000 ml distilled water containing 4 ml glycerol. Shake well a& heat if necessary to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and dispense as desired.

## **Quality Control**

Physical Appearance Cream to light yellow homogeneous free flowing powder

Colour and Clarity of prepared medium Light amber coloured clear solution without any precipitate

**Reaction** Reaction of 4.76% w/v aqueous solution at 25°C. pH : 7.2±0.2

**pH Range** 7.00-7.40





#### **Cultural Response/Characteristics**

DM 2250: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.(Growth evident by turbidity)

Organism	Inoculum (CFU)	Growth	Gas
Escherichia coli ATCC 23724	50-100	good	
Escherichia coli ATCC 39403	50-100	good	
Escherichia coli ATCC 47014	50-100	good	
Escherichia coli ATCC 53868	50-100	good	

#### Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label. **Prepared Media:** 2-8<sup>0</sup> in sealable plastic bags for 2-5 days.

### **Further Reading**

1. Tartoff K. D. and Hobbs C. A., 1987, Improved Media for Growing Plasmid and Cosmid Clones, Bethesda Res. Lab. Focus, 9:12.

 Sambrook J., Fritsch, E. E., and Maniatis T., 1989, Molecular Cloning : A Laboratory Manual, 2nd Ed., Cold Spring Harbor Lab., Cold Spring Harbor, N.Y.

#### **Disclaimer**:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specificatons for identity and performens parameters.

