

## Technical Information

### Tryptone Yeast Extract Agar

**Product Code: DM 2272**

**Application:** - Tryptone Yeast Extract Agar is recommended for estimation of microbial counts in water.

#### Composition\*\*

Ingredients	Gms / Litre
Casein enzymic hydrolysate	6.000
Yeast extract powder	3.000
Agar	12.000
Final pH ( at 25°C)	7.2±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

Periodic sampling and determination of microbial counts of water used for recreation such as beaches etc, open natural/ man made reservoir is important. The total count might be indicative of the overall sanitary conditions pertaining to that site <sup>(1)</sup>. Tryptone Yeast Extract Agar is used as described by ISO Committee <sup>(2)</sup> for the enumeration of viable microorganisms in water under the specification ISO 6222:1988 from these recreation sites

Casein enzymic hydrolysate and yeast extract provide nitrogenous compounds, vitamin B complex and other essential growth nutrients. Total colony forming units (CFU) from the water samples under test is obtained either by spread plate or by pour plate method.

#### Methodology

Suspend 21 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.

#### Quality Control

##### Physical Appearance

Cream to yellow homogeneous free flowing powder

##### Gelling

Firm, comparable with 1.2% Agar gel.

##### Colour and Clarity of prepared medium

Light amber coloured clear to slightly opalescent gel forms in Petri plates.

##### Reaction

Reaction of 2.1% w/v aqueous solution at 25°C. pH : 7.2±0.2

**pH Range:-** 7.00-7.40

##### Cultural Response/Characteristics

DM 2272: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth	Recovery
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	Luxuriant	>=70%
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	>=70%
<i>Salmonella Typhimurium</i> ATCC 14028	50-100	luxuriant	>=70%
<i>Salmonella Enteritidis</i> ATCC 13076	50-100	luxuriant	>=70%



Dehydrated Culture Media  
Bases / Media Supplements

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8<sup>0</sup> in sealable plastic bags for 2-5 days.

## Further Reading

1. Corry J. E. L., Curtis G. D. W., and Baird R. M., Culture Media for Food Microbiology, Vol. 34, Progress in Industrial Microbiology, 1995, Elsevier, Amsterdam
2. International Organization for Standardization (ISO), 1988, Draft ISO/DIS 6222.

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.

