

Technical Information

Dextrose Agar

Product Code: DM 1084

Application: Dextrose Agar is used for cultivation of wide variety of microorganisms and for preparing Dextrose Blood Agar.

Composition**

Ingredients	Gms / Litre
Tryptose	10.000
Beef extract	3.000
Dextrose	10.000
Sodium chloride	5.000
Agar	15.000
Final pH (at 25°C)	7.3±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Dextrose in culture media serves as a source of energy. A basal media with 0.5 - 1.0% dextrose, supplemented with defibrinated blood is recommended for the isolation of a wide variety of fastidious organisms ⁽¹⁾. Dextrose Agar, recommended by APHA ⁽²⁾, contains 1.0% dextrose and therefore supports early and luxuriant growth of a variety of organisms including older cultures. The lag phase is comparatively reduced on this medium. But due to high concentrations of dextrose, the medium is not recommended for studying the haemolytic pattern of organism since dextrose interferes with the haemolytic reaction. Dextrose Agar contains high concentration of dextrose as an energy source for the rapid growth of microorganisms. However this medium is not very suitable for the study of haemolysis because of high carbohydrate content. Beef extract and tryptose serve as sources of nitrogenous compounds, sulphur, carbon, vitamins and minerals. Osmotic balance of the medium is maintained by sodium chloride.

Methodology

Suspend 43 grams of powder media in 1000 ml distilled water. Shake well and heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. If desired, Blood Agar can be prepared by the addition of 5% v/v sterile, defibrinated sheep blood into sterile Dextrose Agar, cooled to 50°C. Mix well and dispense as desired.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Basal medium :Light yellow After addition of 5%v/v sterile defibrinated blood :Cherry red coloured, Basal medium :clear to slightly opalescent gel; After addition :opaque gel forms in Petri plates

Reaction

Reaction of 4.3% w/v aqueous solution at 25°C. pH : 7.3±0.2

pH Range 7.10-7.50

Cultural Response/Characteristics

DM1084: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Growth w/blood	Recovery w/ Blood
<i>Bordetella pertussis</i> ATCC8467	50-100	Good	50-70%	Luxuriant	>=70%
<i>Neisseria meningitidis</i> ATCC 13090	50-100	Good	50-70%	Luxuriant	>=70%
<i>Neisseria gonorrhoeae</i> ATCC 19424	50-100	Good	50-70%	Luxuriant	>=70%
<i>Streptococcus pyogenes</i> ATCC 19615	50-100	Good	50-70%	Luxuriant	>=70%
<i>Clostridium perfringens</i> ATCC 12919	50-100	good	50-70%	luxuriant	>=70%

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8⁰ in sealable plastic bags for 2-5 days.

Further Reading

1. Norton, 1932, J. Lab. Clin. Med., 17:585.
2. Vanderzant C. and Splittstoesser D. F. (Eds.), 1992, Compendium of Methods for the Microbiological Examination of Foods, 3rd Ed., APHA, Washington, D.C.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.