



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Pseudomonas Agar (For Fluorescein)

Product Code: DM 1120

Application: Pseudomonas Agar (For Fluorescein) is recommended for the detection of fluorescein production by *Pseudomonas* species.

Composition**

Ingredients	Gms / Litre
Casein enzymic hydrolysate	10.000
Proteose peptone	10.000
Dipotassium phosphate	1.500
Magnesium sulphate	1.500
Agar	15.000
Final pH (at 25°C)	7.0±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Pseudomonas Agar (For Fluorescein) is based on the formula devised by King et al ⁽¹⁾ and is a modification of U.S. Pharmacopeia ⁽²⁾ for the detection of fluorescein production a water soluble, chloroform insoluble fluorescent pigment by *Pseudomonas* species ⁽³⁾. The medium enhances the elaboration of fluorescein by *Pseudomonas* and inhibits the pyocyanin formation. The fluorescein pigment diffuses from the colonies of *Pseudomonas* into the agar and shows yellow fluorescent colouration. Some *Pseudomonas* strains produce small amounts of pyocyanin resulting in a yellow-green colouration.

Casein enzymic hydrolysate and proteose peptone provide the essential nitrogenous nutrients, carbon, sulphur and trace elements for the growth of *Pseudomonas*. Dipotassium phosphate buffers the medium while magnesium sulphate provides necessary cations for the activation of fluorescein production. Salt concentration exceeding 2% affects pigment production. UV illumination may be bactericidal, so make sure that there is good growth before placing culture under UV light ⁽³⁾.

A pyocyanin-producing *Pseudomonas* strain will usually also produce fluorescein. It must, therefore, be differentiated from other simple fluorescent pseudomonads by other means. Temperature can be a determining factor as most of fluorescent strains fail to grow at 3 5°C. but grow well at 25-30°C ⁽³⁾.

Methodology

Suspend 38 grams of powder media in 1000 ml distilled water containing 10 ml glycerol. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Gelling

Firm, comparable with 1.5% Agar gel

Colour and Clarity of prepared medium

Yellow coloured clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 3.8% w/v aqueous solution (containing 1% v/v glycerol) at 25°C. pH : 7.0±0.2

pH Range 6.80-7.20





Dehydrated Culture Media
Bases / Media Supplements

Cultural Response/Characteristics

DM 1120: Cultural characteristics observed with added 1% glycerol after an incubation at 35-37°C for 18-24 hours.

Organism	Inoculum (CFU)	Growth	Recovery	Colour of colony
<i>Pseudomonas aeruginosa</i> ATCC 17934	50	Luxuriant	>=70%	greenish yellow
<i>Pseudomonas aeruginosa</i> ATCC 27853	50	Luxuriant	>=70%	greenish yellow
<i>Pseudomonas aeruginosa</i> ATCC 9027	50-100	Luxuriant	>=70%	greenish yellow

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8⁰ in sealable plastic bags for 2-5 days.

Further Reading

1. King, Ward and Raney, 1954, J. Lab. Clin. Med., 44: 301.
2. The United States Pharmacopoeia, 2006, USP29/NF24, The United States Pharmacopeial Convention, Rockville, MD.
3. MacFaddin J., 1985, Media for Isolation-Cultivation-Identification-Maintenance of Medical Bacteria, Vol. I, Williams and Wilkins, Baltimore.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents. Do not use the products if it fails to meet specifications for identity and performance parameters.

