

## Technical Information

### Columbia Broth Base

#### Product Code: DM 1145

**Application:** - Columbia Broth Base is used as a general-purpose medium and also for the cultivation of fastidious organisms.

#### Composition\*\*

Ingredients	Gms / Litre
Peptone, special	10.000
Biopeptone	10.000
Heart infusion powder	3.000
L-Cystine hydrochloride	0.100
Dextrose	2.500
Sodium chloride	5.000
Magnesium sulphate	0.100
Ferrous sulphate	0.020
Sodium carbonate	0.600
Tris (hydroxymethyl) aminomethane	0.830
Tris (hydroxymethyl) aminomethane HCl	2.860
Final pH ( at 25°C)	7.5±0.2

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

Morello and Ellner in 1969 devised a liquid medium for the recovery of microorganisms from blood cultures <sup>(1)</sup>. This medium was based on Columbia Blood Agar Base previously formulated by Ellner et al <sup>(2)</sup>. While studying they found that Columbia Broth was superior to a commonly used general-purpose broth for faster growth of *Staphylococcus aureus*,

*Escherichia coli*, viridans Streptococci and Enterococcus groups. In the composition the increased concentration of cystine provided improved recovery of both aerobic and anaerobic microorganisms from blood specimens. Columbia Broth Base when supplemented with SPS (Sodium Polyanethol Sulphonate), a polyanionic anticoagulant inhibited complement and lysozyme activity, interferes with phagocytosis and inactivate aminoglycosides <sup>(3)</sup>. The presence of CO<sub>2</sub> is found to be stimulatory for many organisms. It is an excellent blood culture medium <sup>(4)</sup>. It differs from the agar base as cornstarch is omitted to reduce opalescence <sup>(1)</sup> and salts have been included.

Medium contains peptone special, biopeptone and heart infusion to support luxurious growth of the organisms. Dextrose is added as a carbon and energy source. The medium is buffered with tris buffer. The addition of salts was found to be beneficial for the recovery of organisms. Cystine is the reducing agent. Magnesium and iron are added to facilitate organism growth.

Tube media should be inoculated with 1 to 2 drops of the liquid specimen using a sterile pipette. Swab specimens may be inserted into the broth after inoculation of the plated media. Liquid media should be reduced by placing the tubes with caps loosened under anaerobic conditions for 18-24 hours prior to inoculation for anaerobic incubation. Alternatively, it can be reduced immediately before use by boiling with caps loosened and cooling to room temperature with tightened caps, before inoculation. Growth in tubes is indicated by presence of turbidity compared to an uninoculated control. If growth appears, cultures should be subcultured onto appropriate media. Addition of SPS is inhibitory to *Neisseria* species, and thus 1.2% gelatin addition may counteract the inhibitory effect.

#### Methodology

Suspend 35.01 grams of powder media in 1000 ml distilled water. Shake well & heat to dissolve the medium completely. If desired, SPS (Sodium polyanethol sulphonate) may be added in a final concentration of 0.0 1%. Dispense into tubes and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

## Quality Control

### Physical Appearance

Cream to yellow homogeneous free flowing powder

### Colour and Clarity of prepared medium

Light amber coloured, clear to slightly opalescent solution, may have a fine precipitate.

### Reaction

Reaction of 3.5% w/v aqueous solution at 25°C. pH : 7.5±0.2

**pH range** 7.30-7.70

### Cultural Response/ characteristics

DM 1145: Cultural characteristics observed after an incubation 35-37°C for 18-48 hours.

Organism	Inoculum (CFU)	Growth
Clostridium perfringens ATCC 12924	50-100	good-luxuriant
Neisseria meningitidis ATCC 13090	50-100	good-luxuriant
Staphylococcus aureus ATCC 25923	50-100	good-luxuriant
Streptococcus mitis ATCC 9811	50-100	good-luxuriant
Streptococcus pyogenes ATCC 19615	50-100	good-luxuriant

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8° in sealable plastic bags for 2-5 days.

## Further Reading

1. Morello J. A. and Ellner P. D., 1969, Appl. Microbiol. 17:68.
  2. Ellner P. D., Stoessel C. J., Drakeford E. and Vasi F., 1966, Am. J. Clin. Pathol., 45:502
  3. Reller, Murray and MacLowry, 1982, Cumitech 1A, Blood cultures II, Coord. Ed., ASM, Washington D.C.
- Isenberg, (Ed.), 1992, Clinical Microbiology Procedures Handbook, Vo I. American Society for Microbiology, Washington, D.C

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.