

Technical Information

Phenol Red Lactose Broth

Product Code: DM 1275I

Application: Phenol Red Lactose Broth is recommended for lactose fermentation studies.

Composition**

Ingredients	Gms / Litre
Peptic digest of animal tissue	10.000
Sodium chloride	5.000
Lactose	10.000
Phenol red	0.018
Final pH (at 25°C)	7.5±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Lactose fermentation is a characteristic feature that plays an important role in the primary isolation of pathogenic enterobacteria. This single test makes possible an immediate presumptive distinction between the intestinal commensals & pathogens of various genera which ferment lactose and those which do not ferment lactose⁽¹⁾. Phenol Red Lactose Broth is formulated as recommended by ISO Committee⁽²⁾ for studying lactose fermentation property by coliforms, which is an important differentiating feature for the members of *Enterobacteriaceae*⁽³⁾. Peptic digest of animal tissue provides nitrogenous compounds and other essential growth nutrients. Phenol red is the pH indicator, which turns yellow in acidic condition. 18-24 hours old pure culture is inoculated and incubated at 35 to 37°C for 18-24 hours or upto 30 hours.

Methodology

Suspend 25.02 grams of powder media in 1000 ml distilled water. Mix well and distribute in tubes with inverted Durhams tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Quality Control

Physical Appearance

Light yellow to pink homogeneous free flowing powder

Colour and Clarity of prepared medium

Red coloured clear solution without any precipitate

Reaction

Reaction of 2.5% w/v aqueous solution at 25°C. pH : 7.5±0.2 pH : 7.5±0.2

pH Range:- 7.30-7.70

Cultural Response/Characteristics

DM 1275I: Cultural characteristics observed after an incubation at 35-37°C for 18-24 hours (longer if necessary)

Organism	Inoculum (CFU)	Growth	Acid	Gas
<i>Escherichia coli</i> ATCC 25922	50-100	luxuriant	Positive reaction, Yellow colour	positive reaction
<i>Enterobacter aerogenes</i> ATCC 13048	50-100	luxuriant	Positive reaction, yellow colour	positive reaction
<i>Klebsiella pneumoniae</i> ATCC 13883	50-100	luxuriant	positive reaction, Yellow colour	positive reaction
<i>Proteus vulgaris</i> ATCC 13315	50-100	luxuriant	negative reaction	negative reaction
<i>Salmonella Typhimurium</i> ATCC 14028	50-100	luxuriant	negative reaction	negative reaction
<i>Shigella flexneri</i> ATCC 12022	50-100	Luxuriant	Negative reaction	Negative reaction



Dehydrated Culture Media
Bases / Media Supplements

Storage and Shelf Life-

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8⁰ in sealable plastic bags for 2-5 days.

Further Reading

1. Collee J. G., Fraser A. G., Marmion B. P., Simmons A., (Eds.), 1996, Mackie and McCartney, Practical Medical Microbiology, 14th Edition, Churchill Livingstone
2. International Organization for Standardization, 1990, Draft ISO/DIS 9308-1.
3. Finegold S. M. and Baron E. J., 1986, Bailey and Scotts Diagnostic Microbiology, 7th Ed., The C.V. Mosby Co., St. Louis.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.

