

Technical Information

Ellners Broth

Product Code: DM 1466

Application: - Ellners Broth is used to induce spore formation in *Clostridium perfringens*.

Composition**

Ingredients	Gms / Litre
Proteose peptone	10.000
Yeast extract	3.000
Starch	3.000
Magnesium sulphate	0.100
Monopotassium phosphate	1.500
Disodium phosphate	50.000
Final pH (at 25°C)	7.8±0.2

**Formula adjusted, standardized to suit performance parameters

Principle & Interpretation

Spores are rarely seen in culture (a diagnostic feature) but can be obtained on Ellners Medium ⁽¹⁾. In practice, the routine characterization of clostridia to species level involves morphological examinations, biochemical tests and identification of specific toxins. All clostridia produce spores but they vary markedly in their readiness to do so. Some of which may require prolonged incubation.

C. perfringens are gram-positive rods, often capsulated. In sugar-containing media, the *Clostridium* rods are shorter whereas in protein-containing media, they may become filamentous. Spores formed are usually in small numbers and are not formed in the presence of fermentable carbohydrates. Typically oval, sub-terminal or central spores are formed and are not bulging. Special media like Ellners Broth are used to produce spores ^(1, 2).

Medium is composed of proteose peptone and yeast extract, which supply the necessary nutrients for the growth of the *Clostridia*. Generally sporulation is stimulated by an carbohydrate source and hence starch is included in the medium. Sulphate and phosphate not only buffer the medium but also help in sporulation. *Clostridia* are anaerobic organisms and hence anaerobiosis may be ensured by heating the medium at 100°C for 10 minutes and cooling just before inoculation. It is important that the inoculum should be adequate. 0.5 ml of an actively growing 4-12 hours Meat Broth culture should be introduced with a pipette into the bottom of the tubed medium and incubated anaerobically.

Methodology

Suspend 67.6 grams of powder media in 1000 ml distilled water. Shake well & heat if necessary to dissolve the medium completely. Dispense in tubes. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Quality Control

Physical Appearance

Cream to yellow homogeneous free flowing powder

Colour and Clarity of prepared medium

Amber coloured, clear to slightly opalescent solution

Reaction

Reaction of 6.76% w/v aqueous solution at 25°C. pH : 7.8±0.2

pH range 7.60-8.00



Dehydrated Culture Media
Bases / Media Supplements

Cultural Response/ characteristics

DM 1466: Cultural characteristics observed under anaerobic condition, after an incubation at 35-37°C for 24-76 hours.

Organism	Inoculum (CFU)	Growth	Sporulation
Clostridium perfringens ATCC 12924	50-100	luxuriant	positive

Storage and Shelf Life

Dried Media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

Prepared Media: 2-8° in sealable plastic bags for 2-5 days.

Further Reading

1. Collee J. G., Duguid J. P., Fraser A. G., Marmion B. P., (Eds.), Mackie and McCartney, Practical Medical Microbiology, 1989, 13th Edition, Churchill Livingstone.
2. Trevor W. A., 1977, Anaerobic Bacteriology, 3rd Ed., Butterworths and Co. Ltd.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Donot use the products if it fails to meet specificatons for identity and performens parameters.

