

## Technical Information

### Park and Sanders Broth Base

#### Product Code: DM 2185I

**Application:** - Park and Sanders Broth Base is recommended for selective enumeration of thermotolerant *Campylobacter* species from foods.

#### Composition\*\*

| Ingredients                    | Gms / Litre |
|--------------------------------|-------------|
| Casein enzymic hydrolysate     | 10.000      |
| Peptic digest of animal tissue | 10.000      |
| Yeast extract                  | 2.000       |
| Dextrose                       | 1.000       |
| Sodium citrate                 | 1.000       |
| Sodium chloride                | 5.000       |
| Monohydrogen sodium sulphite   | 0.100       |
| Sodium pyruvate                | 0.250       |
| Final pH ( at 25°C)            | 7.0±0.2     |

\*\*Formula adjusted, standardized to suit performance parameters

#### Principle & Interpretation

Park and Sanders Broth Base (DM2185I) is a slight modification of Park and Sanders Enrichment Broth (DM2185) recommended by APHA <sup>(2)</sup>, for selective enumeration of thermotolerant *Campylobacter* species in food and animal feed.

Casein enzymic hydrolysate, peptic digest of animal tissue, yeast extract provide nitrogenous compounds, carbon, sulphur, vitamins and trace elements. Glucose is the energy source. *Campylobacter* species are microaerophilic. Sodium pyruvate helps for aerotolerance. Sodium sulphite helps in survival of the organism in higher nitrogen atmosphere <sup>(3)</sup>. Sodium citrate suppresses coliforms and gram-positive bacteria. Supplementation of base with antibacterial and antifungal agents are same as described by Park and Sanders <sup>(1)</sup> which markedly reduced growth of normal enteric bacteria and improved recovery of *Campylobacter* species.

After addition of blood and Sanders Selective Supplement I (MS2104), the medium is incubated at 31 to 32°C for 4 hours for the recovery of injured cells. The resuscitation and enrichment of culture must be performed in a microaerobic environment. The organism is sensitive to oxygen and storage at room temperature. The food sample should be stored in an oxygen-free environment with 0.01 % sodium bisulfite and held under refrigeration. Under these conditions the organism will survive 10 times longer than when the same strain is held in a bisulfite-free medium exposed to air at 25°C <sup>(4)</sup>. Some strains of normal enteric organisms may be encountered that are not inhibited or only partially inhibited on this medium.

#### Methodology

Suspend 29.35 grams of powder media in 940 ml distilled water. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Cool to 45°C and aseptically add 50 ml of sterile defibrinated lysed horse blood and reconstituted contents of 1 vial of Park and Sanders Selective Supplement 1 (MS2104). Mix well. Inoculate with food samples and incubate at 31 to 32°C (to recover injured cells) for 4 hours.

Aseptically add reconstituted contents of 1 vial of Park and Sanders Selective Supplement II (MS2105) and incubate at 37°C for 2 hours, then at 42°C under a microaerobic atmosphere for additional 40 to 42 hours with agitation at 100 rpm.

## Quality Control

### Physical Appearance

Light yellow to beige homogeneous free flowing powder

### Colour and Clarity of prepared medium

Basal medium - Light yellow coloured clear solution. After addition of 5% v/v sterile defibrinated lysed horse blood - Cherry red coloured opalescent solution in tubes

### Reaction

Reaction of 2.94% aqueous solution at 25°C. pH : 7.0±0.2

**pH Range** 6.80-7.20

### Cultural Response/Characteristics

DM 2185I: Cultural characteristics observed with added 5% defibrinated lysed horse blood along with MS2104 and MS2105, after an incubation at 42°C for 48 hours under microaerobic atmosphere.

| Organism                               | Inoculum (CFU) | Growth         |
|--|----------------|----------------|
| <i>Campylobacter coli</i> ATCC 33559   | 50-100         | good           |
| <i>Campylobacter jejuni</i> ATCC 29428 | 50-100         | Good-luxuriant |
| <i>Escherichia coli</i> ATCC 25922     | $\geq 10^3$    | inhibited      |

## Storage and Shelf Life

**Dried Media:** Store below 30°C in tightly closed container and use before expiry date as mentioned on the label.

**Prepared Media:** 2-8<sup>0</sup> in sealable plastic bags for 2-5 days.

## Further Reading

1. Park C.E. and Sanders G.W., 1989, Abstr. 5th International Workshop on Campylobacter Infections, Puerto Vallarta, Mexico.
2. Downes F. P. and Ito K., (Eds.), 2001, Compendium of Methods for the Microbiological Examination of Foods, 4th Ed., APHA, Washington, D.C
3. Koidis P. and Doyle M.P., 1983, Eur. J. Clin. Microbiol., 2:384.

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.