



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

G. Vaginalis Selective Supplement

Product Code: MS2056

Application: An antibiotic supplement recommended for the selective isolation of *Gardnerella vaginalis*.

Composition**

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Gentamicin sulphate	2mg
Nalidixic acid	15mg
Amphotericin B	1mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 500 ml of sterile, molten, cooled (45-50°C) Columbia Blood Agar Base DM1144 / DM1144A / Columbia Blood Agar Base, MiVeg VM1144 / VM1144A along with 5% v/v sterile defibrinated sheep blood. Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.