



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Urea Solution

Product Code: MS2253

Application: Filter sterilized urea solution is recommended for the isolation of *Ureaplasma urealyticum* and other *Ureaplasma* species.

Composition**

Per vial sufficient for 1000 ml medium

Ingredients	Concentration
Urea	0.107g
Distilled water	5ml

Directions

Warm up the refrigerated Urea solution to room temperature and aseptically add 1 vial to 790 ml of sterile, molten, cooled (45-50°C) A-7 Agar Base DM2739 along with 195 ml sterile Horse serum BA3239 and 1 vial each of Ureaplasma Selective Supplement MS2254 and Ureaplasma Growth Supplement MS2255 . Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.