



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

ECO157:H7 Selective Supplement

Product Code: MS2247

Application: An antibiotic supplement recommended for selective isolation of *E. coli* O157:H7.

Composition**

Per vial sufficient for 1000 ml medium

*Ingredients	Concentration
Vancomycin	8mg
Cefixime	0.050mg
Cefsulodin	10mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml sterile distilled water. Mix gently to dissolve the contents completely. Aseptically add the rehydrated contents to 1000 ml of sterile, cooled to 45-50°C Modified Soyabean Bile Broth Base DM1286I / Buffered Peptone Water DM1614 or Buffered MiVeg Peptone Water VM1614 . Mix well and dispense into sterile tubes.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.