



Dehydrated Culture Media  
Bases / Media Supplements

## Technical Information

### Legionella Selective Supplement II

**Product Code: MS2037**

**Application:** An antibiotic supplement recommended for the selective isolation of *Legionella* species from mixed cultures.

#### Composition\*\*

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Cephalothin	2mg
Colistin sulphate	8mg
Vancomycin	0.250mg
Cycloheximide	40mg

#### Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 500 ml of sterile, molten, cooled (45-50°C) Buffered Charcoal Yeast Extract Agar Base DM1813 / Modified Buffered Charcoal Agar Base DM1892 / Modified Buffered Charcoal MiVeg Agar Base VM1892 alongwith rehydrated contents of Legionella Selective Supplement IV MS2040 and Legionella Supplement MS2041. Mix well and pour into sterile petri plates.

#### Storage and Shelf Life

Store at 2 - 8°C. Use before expiry date on the label.

#### Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.