



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Ticarcillin Supplement

Product Code: MS2102

Application: An antibiotic supplement recommended for the selective enrichment of *Yersinia enterocolitica*.

Composition**

Per vial sufficient for 1000 ml medium

*Ingredients	Concentration
Ticarcillin	1mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 988 ml of sterile, molten cooled (45-50°C) ITC Broth Base DM2220 / ITC MiVeg Broth Base VM2220 containing 10 ml Potassium Chlorate Supplement MS2103 . Mix well and dispense into sterile tubes.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.