

## Technical Information

### MiCrome Nickels & Leesment Selective Supplement

**Product Code: MS2245**

**Application:** An antibiotic supplement used for the selective isolation of citrate-fermenting lactic acid bacteria.

#### Composition\*\*

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Vancomycin	100mg

#### Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and aseptically add it to sterile, molten, cooled (45-50°C) MiCrome Nickels and Leesment Medium DM2712 . Mix well and pour into sterile Petri plates.

#### Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

#### Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.