



Dehydrated Culture Media  
Bases / Media Supplements

## Technical Information

### Kimmig Selective Supplement

**Product Code: MS2111**

**Application:** An antimicrobial supplement recommended for selective isolation and cultivation of fungi.

#### Composition\*\*

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Part A	
Colistin sulphate	40mg
Part B	
Novobiocin	50mg

#### Directions

Rehydrate the contents of one vial each of Part A and Part B, separately with 2.5 ml sterile distilled water and mix well. On reconstitution aseptically add to 500 ml sterile, molten, cooled (45-50°C) Kimmig Fungi Agar Base DM2232 / Kimmig Fungi MiVeg Agar Base VM2232 . Mix well and pour into sterile petri plates.

#### Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

#### Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.