



Dehydrated Culture Media  
Bases / Media Supplements

## Technical Information

### Vancomycin Supplement

**Product Code: MS2261**

**Application:** This antibiotic supplement is recommended for selective isolation of Vancomycin Resistant Enterococci.

#### Composition\*\*

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Vancomycin	3mg

#### Directions

Rehydrate the contents of 1 vial with 5 ml of sterile distilled water and aseptically add 2 vials to 1000 ml of sterile molten Vancomycin Resistant Enterococci (VRE) Agar Base DM2763 along with 1 vial of Meropenem Supplement MS2262 .

Mix well and pour into sterile Petri plates.

#### Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

#### Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.