



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Park and Sanders Selective Supplement I

Product Code: MS2104

Application: An antibiotic supplement recommended for the selective isolation and enumeration of thermotolerant *Campylobacter* species

Composition**

Per vial sufficient for 950 ml medium

*Ingredients	Concentration
Vancomycin	10mg
Trimethoprim	10mg

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of 50% ethanol. Mix well and aseptically add to 940 ml of sterile, cooled Park and Sanders Enrichment Broth Base DM2185 / Park and Sanders Broth Base DM2185I / Park and Sanders Enrichment MiVeg Broth Base VM2185 along with 50 ml sterile defibrinated blood and Park and Sanders Selective Supplement II MS2105 . Mix well and dispense into sterile tubes.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.