



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

BSK-H Supplement

Product Code: MS2180

Application: A sterile filtered rabbit serum to support growth of *Borrelia burgdorferi*

Composition**

Ingredients	Concentration
rabbit serum	100ml

Directions

To prepare complete medium warm up the refrigerated rabbit serum to room temperature and aseptically add 30 ml to 500 ml sterile, cooled (45-50°C) BSK-H Medium Base DM2668 / LQ001 / BSK-H Medium Powder DM6107 along with one vial of Antibiotic Mixture of *Borrelia* MS2179 . Mix well and dispense as desired.

Storage and Shelf Life

Store at -20°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.