

## Technical Information

### CAT Selective Supplement

**Product Code: MS2145**

**Application:** An antibiotic supplement recommended for selective isolation of thermophilic *Campylobacter* species.

### Composition\*\*

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Cefoperazone	4mg
Teicoplanin	2mg
Amphotericin B	5mg

### Directions

Rehydrate the contents of one vial aseptically with 2 ml of sterile distilled water. Mix well and add aseptically to 500 ml of sterile, molten, cooled (45-50°C) Blood Free Campylobacter Selectivity Agar Base DM1887 / Blood Free Campylobacter Selectivity MiVeg Agar Base VM1887. Mix well and pour into sterile petri plates.

### Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

### Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.