



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Trichomonas Selective Supplement I

Product Code: MS2099

Application: An antimicrobial supplement recommended for selective isolation of *Trichomonas* species from clinical specimens.

Composition**

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Streptomycin	500mg
Penicillin G	125000Unit

Directions

Rehydrate the contents of 1 vial aseptically with 5 ml of sterile distilled water. Mix well and avoid frothing of solution. Aseptically add the contents to 475 ml of sterile, molten, cooled (45-50°C) Kupferberg Trichomonas Broth Base (Trichomonas Broth Base, Kupferberg) DM1305 / Kupferberg Trichomonas MiVeg Broth Base VM1305 along with 25 ml sterile beef or human serum. Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.