

Technical Information

CNA Supplement

Product Code: MS2115

Application: An antibiotic supplement recommended for selective isolation of anaerobic gram-positive cocci.

Composition**

Per vial sufficient for 950 ml medium

*Ingredients Concentration
Colistin sulphate 10mg
Nalidixic acid 10mg

Directions

Rehydrate the contents of 1 vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add to 950 ml of sterile, molten, cooled (45-50°C) Schaedler Agar DM1291 / Schaedler MiVeg Agar VM291 along with 50 ml sterile defibrinated sheep blood and Vitamin K1 Supplement MS2114. Mix and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at CDH is true and accurate.
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.