

Technical Information

Bifidobacterium Selective Supplement

Product Code: MS2285

Application: An supplement recommended for the selective isolation of *Bifidobacterium* species.

Composition**

Per vial sufficient for 1000 ml medium

*Ingredients	Concentration
Propionic acid	5ml

Directions

Warm up the refrigerated contents of one vial to room temperature. Mix well and aseptically add the contents to 1000ml of sterile, molten, cooled (45-50°C) Bifidobacterium Agar Modified DM2858 . Mix well and pour into sterile Petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.