

Technical Information

IMRV/RV Selective Supplement

Product Code: MS2193

Application: An antibiotic supplement recommended for isolation of *Salmonella* from food stuffs and other materials.

Composition**

Per vial sufficient for 500ml / 1000ml medium

*Ingredients	Concentration
Novobiocin	10mg

Directions

Rehydrate the contents of one vial aseptically with 5 ml of sterile distilled water and aseptically add it to 1000 ml sterile, molten, cooled (45-50°C) Semisolid IMRV Medium Base DM2427 / Semisolid IMRV MiVeg Medium Base VM2427 / Modified Semisolid RV Medium Base DM2482 & 500 ml of Semisolid RV Medium Base DM2428 / Semisolid RV MiVeg Medium Base VM2428 . Mix well and pour into sterile petri plates.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.