



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

Oxacillin Resistance Selective Supplement

Product Code: MS2191

Application: An antibiotic supplement recommended for screening oxacillin resistant microorganisms.

Composition**

Per vial sufficient for 500 ml medium

*Ingredients	Concentration
Oxacillin	1mg
Polymyxin B	25000IU

Directions

Rehydrate the contents of one vial aseptically with 2 ml of sterile distilled water. Mix well and aseptically add it to 500 ml sterile, molten, cooled (45-50°C) Oxacillin Resistance Screening Agar Base DM2454 . Mix well and dispense as desired.

Storage and Shelf Life

Store at 2-8°C. Use before the expiry date on the label.

Disclaimer

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.