

## Technical Information

### CNV Supplement

#### Product Code: PHS1014

**Application:** Selective supplement for the detection of *Pseudomonas savastanoi pv. phaseolicola*, *Pseudomonas syringae* and *Xanthomonas axonopodis* in seeds of beans.

#### Formula

Per vial sufficient 1000ml medium

Cephalexin monohydrate	80mg
Nystatin	35mg
Vancomycin Hydrochloride	10mg

#### Directions

Rehydrate the contents of 1 vial with 2 ml of 0.2N NaOH. Mix well and aseptically add to sterile molten Phyto Tyrosine Peptone Agar Base (MT Medium) (PHM1018). Mix well and pour into sterile Petri plates.

#### Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performs parameters.