

## Technical Information

### Phyto Lactose Trehalose Agar Base

**Product Code: PHM1006**

**Application:** Semi-selective medium for the detection *Xanthomonas hortorum pv carotae* in carrot.

### Composition\*\*

Ingredients	Grams/Litre
Yeast extract	0.50
Potassium dihydrogen phosphate	1.2
Di-potassium hydrogen phosphate	1.2
Ammonium chloride	1.00
Lactose monohydrate	10.00
Trehalose monohydrate	4.00
2-Thiobarbituric acid	0.20
Agar	18.00
<b>Total</b>	<b>36.10 gm/liter</b>
Final pH (at 25°C)	6.6

\*\*Formula adjusted standard to suit the performance parameter

### Principle And Interpretation

The causative agent of bacterial leaf blight of carrot is *Xanthomonas hortorum pv. carotae* (also known as *Xanthomonas campestris pv. carotae*) causes significant yield losses. For routine testing of carrot seed a combination of two semi-synthetic media, Phyto Lactose Trehalose Agar Base/ Phyto Xanto Carota Agar Base and Phyto Lactose Trehalose Agar Base/ Phyto peptone Bromide Agar Base .is recommended (1). This medium is a modification of the KM-1 medium (2). This medium is more sensitive and has a greater reproducibility as compared to other media.

Yeast extract serves as a nutritional source. Phosphates buffers the medium. Nystatin in the medium is used to inhibit the growth of fungi. Lactose monohydrate and trehalose monohydrate are the fermentable carbohydrates which are an additional nutritional source. This medium contains lowered concentration of tobramycin (2mg) since some strains of *Xanthomonas hortorum pv. carotae* are sensitive to concentration of tobramycin sulphate greater than 4 mg.

### Directions

Suspend 36.10 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Sterilize the medium by autoclaving at 15 lbs pressure (121°C) for 15 minutes .Cool to 45-50°C and aseptically add the rehydrated contents of one vial of CNTB supplement-1 (PHS1005). Mix well and pour into sterile Petri plates .

### Quality Control

**Appearance:**

Cream to yellow coloured, homogeneous, free flowing powder.

**Gelling:**

Firm, comparable with 1.8% Agar gel.

**Colour and Clarity of prepared medium:**

Yellow coloured, opalescent gel with white precipitate forms in Petri plates

**Reaction:**

Reaction of 3.51% w/v aqueous solution is pH 6.6 at 25°C.

**Cultural Response:**

Cultural characteristics observed with added CNTB supplement-1 (PHS1005), after incubation at 30-32°C for 5-6 days.

### Organism (ATCC)

### Growth

### Colony characteristics

*Xanthomonas hortorum pv carotae*

Luxuriant

Cream to yellow or light brown or peach yellow, glistening, Round and about 2-4 mm in diameter.

*Staphylococcus aureus* (25923)

Inhibited

-

## Storage and Shelf Life

Store below 30°C and the prepared medium at 2 - 8°C. Use before expiry date on the label.

## Further Reading

- Detection of *Xanthomonas hortorum pv. carotae* on *Daucus carota*. International Rules for seed testing 7-020, 2006.
- Kim, H.K., Sasser, M. and Sands, D.C. 1982. Selective medium for *Xanthomonas hortorum pv translucens*. Phytopathology 72:936.

## Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.