



Dehydrated Culture Media
Bases / Media Supplements

Technical Information

BME Vitamin Solution 100X Sterile filtered

Product Code : VA1001

Application:-

VA1001 is a 100X concentrate of various vitamins present in Basal Medium Eagle (BME). Addition of vitamin supplement to medium stimulates growth and prolongs viability of cells in culture.

VA1001 is a sterile filtered BME vitamin solution prepared in 0.85% normal saline.

Composition**

Ingredients	mg / Litre
VITAMINS	
Choline chloride	100.000
D-Biotin	100.000
D-Ca-Pantothenate hemicalcium	100.000
Folic acid	100.000
Niacinamide	100.000
Pyridoxal hydrochloride	100.000
Riboflavin	10.000
Thiamine hydrochloride	100.000
i-Inositol	200.000

Methodology

For preparing 1X concentration of vitamin solution use 10ml of 100X solution of VA1001 per litre of medium i.e. 10ml of 100X solution VA1001 contains 1X vitamins.

Quality Control

Appearance

Yellow colored, clear solution.

pH

7.20 -7.80

Osmolality in mOsm/Kg H₂O

270.00-310.00

Sterility

No bacterial or fungal growth is observed after 14 days of incubation, as per USP specification.

Cultural Response

The growth promotion capacity of the solution is assessed qualitatively by analyzing the cells for the morphology and quantitatively by estimating the cell counts and comparing it with a control solution.

Endotoxin content

NMT 5EU/ml



Dehydrated Culture Media

Storage and Shelf Life

Store at -20°C.

Shelf life of the product is 12 months.

Use before expiry date given on the product label.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.