

Technical Information

M9 Minimal Medium Salts (5X)

Product Code: G1013

M9 Minimal Medium Salts (5X) is recommended for growing *Escherichia coli* for molecular biology purpose

Composition**

Ingredients	Grams/Litre
Disodium hydrogen phosphate	33.90
Potassium dihydrogen phosphate	15.00
Sodium chloride	2.50
Ammonium chloride	5.00

** Formula adjusted, standardized to suit performance parameters

Methodology

Suspend 56.4 grams in 1000 ml distilled water. Heat to boiling to dissolve the medium completely. Dispense as desired and sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes.

Principle and Interpretation

M9 Minimal Salts, 5x is used in preparing M9 Minimal medium used for the cultivation and maintenance of *E. coli* in molecular biology. M9 is a chemically defined minimal growth medium for *E. coli* cultures (1, 2). It can be supplemented with specific amino acids or other required nutrients for the selection of some specific auxotroph. Sometimes the incorporation of certain additives (e.g. thiamine or casamino acids) enhances the bacterial growth. In 1959, Jacob and Monod published an article where they discussed about minimal media of *E. coli* (3). The components of this media supply the required nutrients for proper growth. Ammonium chloride acts as the nitrogen source and Sodium chloride maintains the final osmolarity of the medium. M9 Minimal Salts Base, 5x, is a 5x concentrate and should be diluted to 1x concentration before usage. Glucose may be added as a carbon source. Some other nutritional elements (e.g. calcium and magnesium) can be added as supplement.

Quality control

Appearance of Powder :

White to cream coloured, homogeneous, free flowing powder.

Colour and Clarity :

Colourless, clear solution without any precipitate.

Cultural Response :

Cultural characteristics observed after an incubation at 35-37°C for 18 - 48 hours.

Organisms (ATCC)

Escherichia coli ATCC 23724

Escherichia coli ATCC 25922

Escherichia coli MTCC 1652

Growth

good-luxuriant

good-luxuriant

good-luxuriant

Storage and Shelf Life

Store below 30°C and the prepared medium at 2 - 8°C. Use before expiry date on the label.

Reference

1. Davis, L.G., M.D. Dibner and J.F. Battey, Basic methods in molecular biology, Elsevier, new York, (1986).
2. Sambrook, J., E. F. Fritsch, and T. Maniatis, 1989, Molecular cloning: a laboratory manual, 2nd edition ed., Cold Spring Harbour laboratory, Cold Spring Harbour, N.Y.
3. Pardee, A. B., F. Jacob, and J. Monod. 1959. The genetic control and cytoplasmic expression of "inducibility" in the synthesis of β -galactosidase in *E. coli*. J. Mol. Biol. 1:165-178.

Disclaimer :

- User must ensure suitability of the product(s) in their application prior to use.
- The product conforms solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate.
- **Central Drug House Pvt. Ltd.** reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specifications for identity and performance parameters.