

Bases / Media Supplements

Technical Information

ISP Medium No. 5 (Glycerol Asparagine Agar Base)

Product Code: DM 1360

Application: - ISP Medium No. 5 (Glycerol Asparagine Agar Base) is recommended for cultivation of *Streptomyces* species as per International Streptomyces Project.

Composition**		
Ingredients	Gms / Litre	
L-Asparagine	1.000	
Dipotassium phosphate *Trace salt solution (ml) Agar 1ml of Trace salt solution contains	1.000 1.000 20.000	
Ferric sulphate heptahydrate	0.001	
Manganese chloride tetrahydrate	0.001	
Zinc sulphate heptahydrate	0.001	
Final pH (at 25°C) **Formula adjusted, standardized to suit perfor parameters	7.4±0.2 mance	

Principle & Interpretation

ISP Medium No. 5 (Glycerol Asparagine Agar Base) has composition described by Shirling and Gottlieb ⁽¹⁾ and is used for cultivation and characterization of *Streptomyces* species as recommended by the International Streptomyces Project. Being primarily soil inhabitants, *Streptomyces* are one of the most common agents of causing actinomycotic mycetoma. Areas of body more prone to formation of mycetomas are those that are frequently traumatized or come in contact with soil ⁽²⁾.

This medium provides consistent and reproducible characteristic features which help in the identification of Streptomyces. Glycerol serves as the carbon source while asparagine is the amino acid source for the growth of Streptomyces species. Trace mineral requirement of Streptomyces is satisfied by the trace salt solution, which contains various salts. Dipotassium phosphate buffers the medium.

Methodology

Suspend 23 grams of powder media in 1000 ml distilled water containing 10 ml glycerol. Shake well & heat to dissolve the medium completely. Sterilize by autoclaving at 15 lbs pressure (121°C) for 15 minutes. Mix well and pour into sterile Petri plates.

Quality Control

Physical Appearance

Off-white to yellow homogeneous free flowing powder **Gelling** Firm, comparable with 2.0% agar gel.

Colour and Clarity of prepared medium Light amber coloured, clear to slightly opalescent gel forms in Petri plates

Reaction

Reaction of 2.3% w/v aqueous solution containing 1.0% glycerol at 25°C. pH : 7.4 \pm 0.2





Dehydrated Culture Media Bases / Media Supplements

pH Range 7.20-7.60

Cultural Response/Characteristics

DM 1360: Cultural characteristics observed after an incubation at 25-30°C upto 15 days.

Organism		Growth	
Streptomyces albus subsp albus ATCC 3006		good-luxuriant	
Streptomyces lavendulae ATC 8664	С	good-luxuriant	
Streptomyces peucetius ATCC 29050		good-luxuriant	

Storage and Shelf Life

Dried media: Store below 30°C in tightly closed container and use before expiry date as mentioned on the label. **Prepared Media:** 2-8⁰ in sealable plastic bags for 2-5 days.

Further Reading

1. Shirling E. B. and Gottlieb D., 1966, International J. Systemic Bacteriol., 16:3.

2. Murray P. R., Baron J. H., Pfaller M. A., Jorgensen J. H. and Yolken R. H., (Eds.), 2003, Manual of Clinical Microbiology, 8th Ed., American Society for Microbiology, Washington, D.C.

Disclaimer:

- User must ensure suitability of the product(s) in their application prior to use.
- The product conform solely to the technical information provided in this booklet and to the best of knowledge research and development work carried at **CDH** is true and accurate
- Central Drug House Pvt. Ltd. reserves the right to make changes to specifications and information related to the products at any time.
- Products are not intended for human or animal diagnostic or therapeutic use but for laboratory, research or further manufacturing of diagnostic reagents extra.
- Statements contained herein should not be considered as a warranty of any kind, expressed or implied, and no liability is accepted for infringement of any patents.
- Do not use the products if it fails to meet specificatons for identity and performens parameters.

